

MOBILE SURGI-STRETCHER



ELECTRIC UNITS FOR OPHTHALMOLOGY

ELECTRIC EYE-STRETCHER (5E82EYEST)

Designed as a mobile surgical platform for EYE/ENT and plastic surgery procedures.

Key Features:

- Low-profile headpiece allowing site access for both superior (over brow) and lateral (side approach) procedures.
- Battery-Powered with interactive hand pendant.
- Available controls: Height adjustment, Backrest & Leg Section, Trendelenurg/Reverse Trendelenburg.
- Stainless steel Retracto® side rails with rail caps.
- 3-stage Telescoping IV Pole (permanently mounted at the foot-end).
- Emergency Backrest release to override electric/battery.
- Base cover with built-in O₂ tank holder and storage space.



Patented Articulating Headrest
Dual operating knobs for infinite positioning and adjustability.



Ergonomic Push Handles
For ease of handling and safety for both staff and patients.



363
KG



Intuitive Hand Pendant

Easy commands for movement and presets with memory for up to 3 user-created configurations.



Smart Castor Technology

Automatically self-breaks after a short duration of no movement.

Specifications:	5E82EYEST	5E82EYXST
	standard width	wide width
	cm †	cm †
Overall Width:	81	94
Overall Length:	235	235
Overall Height: minimum	56	56
Overall Height: maximum	102	102
Patient Surface: width	64	76
Patient Surface: length	217	217
Side Rail Length:	124	124
Side Rail Height:	36	36
Mattress Pad: thickness	8	9
Casters:	15	15
	degrees [°]	degrees [°]
Backrest: movement	0 - 85	0 - 85
Knee Flex: movement	0 - 60	0 - 60
Trendelenburg: standard	15	15
Trendelenburg: reverse	15	15
	power	power
Volts (V)	100 - 240	100 - 240
Hertz (Hz)	~ 50 / 60	~ 50 / 60
	kg †	kg †
Weight Capacity: EVENLY DISTRIBUTED	363	363
Shipping Weight:	191	209