

## MODEL : EC-3FN

### EXAMINATION ELECTRIC COUCH 3-SECTION



<b>OVERALL DIMENSION</b>	: L2020x W625 x H475~995 mm
<b>MATERIAL</b>	: MILD STEEL FRAMEWORK
<b>FINISHING</b>	: EPOXY POWDER COATING (RAL 9001 BEIGE)
<b>UPHOLSTERY</b>	: PVC / DARK BLUE, MAROON, DARK GREY, BLACK
<b>CASTOR DIAMETER</b>	: Ø5" / 125mm WHEELS WITH LOCKS
<b>POWER SUPPLY VOLTAGE</b>	: 230 VAC / 50 Hz @ 2.0A MAX.
<b>MAX. LOAD CAPACITY</b>	: 220 KG

#### STANDARD FEATURES

- Adjustable Backrest (-30° to 70°), Footrest (0° to -90°) and Hi-Lo function 475mm to 995mm
- Electric actuators by foot pedal and handset control
- 65mm thick foam mattress with PVC covering
- Paper roll holder for up to 60 cm width paper roll

#### OPTIONAL

- Cut out face hole on backrest
- Side rails at middle section on both sides (as shown above)
- Customizable overall dimension
- Customizable backrest, footrest and Hi-Lo adjustment