

MODEL : EC6200-700-ECHO
ELECTRICAL ECHO EXAMINATION COUCH



OVERALL DIMENSION	: L1900 x W700 x H500-850mm
MATERIAL	: MILDSTEEL FRAMEWORK, STAINLESS STEEL RAIL
FINISHING	: EPOXY POWDER COATING (LIGHT GREY TEXTURE)
ECHO CUT-OUT AREA	: L300 x W220 mm (LEFT SIDE)
RETRACTABLE SEAT	: L300 x W250 mm (RIGHT SIDE)
CASTOR DIAMETER	: 125mm SWIVEL CASTOR WITH INDIVIDUAL BRAKE
UPHOLSTERY COLORS AVAILABLE	: BEIGE / YELLOW / ORANGE / RED / MAROON DARK GREEN / BLUE / DARK BLUE / BLACK

STANDARD FEATURES

- HEIGHT ADJUSTMENT 500MM TO 800MM ABOVE FLOOR LEVEL
- BACKREST INCLINATION 0 to 75 DEGRESS
- THICK FOAM MATTRESS 65MM, PVC COVER
- HANDSET CONTROLLER FOR HI-LO & BACKREST AJDUSTMENT
- FOOT CONTROL PEDAL FOR HI-LO & BACKREST AJDUSTMENT
- PAPER ROLL HOLDER (L600mm) STAINLESS STEEL AT HEAD END OR FOOT END
- LEFT SIDE ONLY STAINLESS STEEL SIDERAIL
- GAP BETWEEN 2 SECTIONS LESS THAN10MM
- ECHO CUT-OUT AREA 15CM FROM TOP OF BACKREST SECTION
- RETRACTABLE SEAT 30CM FROM TOP OF BACKREST SECTION
- HEIGHT ADJUSTABLE IV POLE BRACKET AT TOP LEFT CORNER