

MODEL : EC 6200
EXAMINATION ELECTRIC COUCH 2-SECTION



Optional with Side Rails



Optional with Arm Rests

OVERALL DIMENSION	: L2020 x W625 x H590~1035 mm
MATERIAL	: MILDSTEEL FRAMEWORK
FINISHING	: EPOXY POWDER COATING
UPHOLSTER	: PVC / DARK BLUE, MAROON, DARK GREY, BLACK
CASTOR DIAMETER	: ø3" / 75mm SWIVEL CASTOR
POWER SUPPLY VOLTAGE	: 230 VAC / 50 Hz @ 2.0A MAX.
MAX. LOAD CAPACITY	: 300 KG

STANDARD FEATURES

- Adjustable Backrest (-30° to +45°) and Hi-Lo function 590mm to 1035mm
- Electric actuators by foot pedal and handset control
- Stabilize frame through pedal mechanism
- 65mm thick foam mattress PVC covering
- Paper roll holder for up to 60 cm width paper roll

OPTIONAL

- Cut out face hole on backrest
- Side rails at body section on both sides
- Customizable overall dimension
- Customizable backrest, footrest and Hi-Lo adjustment