

INTRODUCING THE NEW
HSS SERIES SURGICAL STOOL



HSS500-S SADDLE SEAT STYLE



HSS500-WD WATER DROP SEAT STYLE



HSS500-W WEDGE SEAT STYLE

BUILT FOR PERFORMANCE
AND COMFORT

HAUSTED SURGICAL STOOL

IDEAL FOR EYE, ENT AND PLASTIC PROCEDURES



Adjustable Contoured Backrest

360° Rotating Urethane Arm Rest for Comfortable Arm Positioning

300 POUNDS
300 lb (136 kg)
Weight Capacity
EVENLY DISTRIBUTED

Dual-density Foam Seat with **DURA-GEL®** Insert for Enhanced Comfort



180° Arm Rest Support for more Stable Arm Positioning

Adjustable Height Range
21 - 28 inches (53 - 71 cm)

OPTIONAL ACCESSORY:
Riser Kit Increases Height Range
an Additional 2 inches (5 cm)
n No: HSS-Riser

ABS Thermoformed Hood for Easy Cleaning



Front Locking System

FOOT-CONTROLLED HYDRAULICS
with Three-Spindle Base and Four 3 inch (8 cm) Dual Tente® Casters for Greater Stability



Complete the HSS Series Surgical Stool Experience with the new Hausted **ESC² Ophthalmology Platform**
Item No: ESC2EYEST

ENGINEERED TO PROVIDE
ENHANCED
SURGEON COMFORT

WARNINGS: Cancer and Reproductive Harm -
www.P66Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



EC REP

CE Partner4U BV
Esdoornlaan 13
3951 DB Maarn
The Netherlands
www.cepartner4u.eu

GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company
MDSAP
MEDICAL DEVICE SINGLE AUDIT PROGRAM
Intertek



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com

graham field

© 2020, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE Casters, Inc.