

INSTRUMENT CLEANING AND MAINTENANCE

The Key Surgical Instrument Cleaning line is designed to provide you with the tools needed to clean your surgical instruments. In addition to our wide variety of cleaning brushes, we also offer enzymatic solutions to help improve the decontamination process.

With bristles available in nylon, brass, polypropylene or stainless steel, we also provide an incredible variety of lengths and diameters to fit your needs.

Key Surgical Endoscope Brushes are available in a wide range of diameters and lengths designed for manually cleaning the channels of a variety of endoscopes such as colonoscopes, gastroscopes, ureteroscopes, and duodenoscopes. Our selection of brushes to clean valves, suction wells, or elevator mechanisms as well as brush kits will have your department set up for success.

With our vast assortment of brushes - from channel cleaning to endoscope cleaning to steriliser cleaning - we can help you match the right brush to the right cleaning task at hand.

Channel Cleaning Brushes

Pipe Cleaners

Toothbrush-Style Cleaning Brushes

General Instrument Cleaning Brushes

Steriliser Cleaning Brushes


Brush Assortments

Endoscope Cleaning Brushes

Please note: Unless specified, all products are not made with natural rubber latex.



Instrument Cleaning – Channel Cleaning Brushes




BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-06-039	1mm	6.35mm	152mm	2 pcs
B	BR-08-039	1mm	12.7mm	203.2mm	2 pcs
	BR-18-039	1mm	6.35mm	457mm	2 pcs
C	BR-06-062	1.58mm	12.7mm	152mm	2 pcs
	BR-18-062	1.58mm	12.7mm	457mm	2 pcs
D	BR-08-070	1.77mm	12.7mm	203.2mm	2 pcs
	BR-18-070	1.77mm	12.7mm	457mm	2 pcs
	BR-24-070	1.77mm	12.7mm	609.6mm	2 pcs
E	BR-06-079	2mm	12.7mm	152mm	2 pcs
F	BR-08-079	2mm	19.05mm	203.2mm	2 pcs
	BR-12-079	2mm	19.05mm	305mm	2 pcs
G	BR-24-080	2mm	25.4mm	609.6mm	2 pcs
H	BR-08-087	2.2mm	12.7mm	203.2mm	2 pcs
	BR-18-087	2.2mm	12.7mm	457mm	2 pcs
I	BR-08-092	2.3mm	12.7mm	203.2mm	2 pcs
J	BR-24-098	2.49mm	12.7mm	609.6mm	2 pcs
K	BR-08-105	2.67mm	25.4mm	203.2mm	2 pcs
	BR-18-105	2.67mm	25.4mm	457mm	2 pcs
L	BR-08-118	3mm	38.1mm	203.2mm	2 pcs
M	BR-12-118	3mm	69.85mm	305mm	2 pcs

Instrument Cleaning – Channel Cleaning Brushes




BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-18-125	3.18mm	50.8mm	457mm	2 pcs
B	BR-08-131	3.3mm	19.05mm	203.2mm	2 pcs
C	BR-08-144	3.66mm	19.05mm	203.2mm	2 pcs
D	BR-08-158	4mm	19.05mm	203.2mm	2 pcs
	BR-12-158	4mm	19.05mm	305mm	2 pcs
	BR-18-158	4mm	19.05mm	457mm	2 pcs
E	BR-12-197	5mm	76.2mm	305mm	2 pcs
	BR-16-197	5mm	76.2mm	406mm	2 pcs
	BR-18-197	5mm	76.2mm	457mm	2 pcs
	BR-24-197	5mm	76.2mm	609.6mm	2 pcs
F	BR-24-200	5.08mm	50.8mm	609.6mm	2 pcs
G	BR-12-236	6mm	50.8mm	305mm	2 pcs
	BR-16-236	6mm	50.8mm	406mm	2 pcs
	BR-24-236	6mm	50.8mm	609.6mm	2 pcs

Instrument Cleaning – Channel Cleaning Brushes



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-06-250	6.35mm	50.8mm	152mm	2 pcs
	BR-08-250	6.35mm	50.8mm	203.2mm	2 pcs
	BR-18-250	6.35mm	50.8mm	457mm	2 pcs
B	BR-16-276	7mm	76.2mm	406mm	2 pcs
	BR-24-276	7mm	76.2mm	609.6mm	2 pcs
C	BR-24-312	7.9mm	76.2mm	609.6mm	2 pcs
D	BR-06-375	9.5mm	38.1mm	152mm	2 pcs
E	BR-24-375	9.5mm	76.2mm	609.6mm	2 pcs
F	BR-16-393	10mm	76.2mm	406mm	2 pcs
	BR-18-393	10mm	76.2mm	457mm	2 pcs
	BR-24-393	10mm	76.2mm	609.6mm	2 pcs
G	BR-24-425	10.8mm	76.2mm	609.6mm	2 pcs

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes



A



B



C



D



E

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-08-440	11.18mm	76.2mm	203.2mm	2 pcs
	BR-24-440	11.18mm	76.2mm	609.6mm	2 pcs
B	BR-16-443	11.25mm	76.2mm	406mm	2 pcs
	BR-24-443	11.25mm	76.2mm	609.6mm	2 pcs
C	BR-16-472	12mm	76.2mm	406mm	2 pcs
D	BR-06-500	12.7mm	50.8mm	152mm	2 pcs
E	BR-16-591	15mm	50.8mm	406mm	2 pcs
	BR-18-591	15mm	50.8mm	457mm	2 pcs

K

Instrument Cleaning

Instrument Cleaning – Channel Cleaning Brushes



A



B



C



D

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-08-630	16mm	50.8mm	203mm	2 pcs
B	BR-16-787	20mm	50.8mm	406mm	2 pcs
C	BR-18-880	22.35mm	88.9mm	457mm	2 pcs
D	BR-24-880	22.35mm	107.95mm	609.6mm	2 pcs

K

Instrument Cleaning

Instrument Cleaning – Channel Cleaning Brushes

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end (pull through style excluded).

Product code	Description	Diameter	Bristle area length	Length	Qty
A BR-3044	nylon	0.5mm	6.35mm	76mm	10 pcs
B BR-3045	nylon	0.5mm	6.35mm	203.2mm	10 pcs
C BR-3047	nylon	0.9mm	12.7mm	101.6mm	10 pcs
D BR-4120	nylon, centered (pull through)	6.35mm to 15.24mm	50.8mm	457.2mm	2 pcs
E BR-4140	nylon, centered (pull through)	12.7mm	76.2mm	609.6mm	6 pcs



A



B




C



D

(PULL THROUGH)



E


(PULL THROUGH)

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

K

Instrument Cleaning

Instrument Cleaning – Channel Cleaning Brushes



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities also available in 50 pcs.

	Product code	Description	Diameter	Bristle area length	Length	Qty
A	BR-8110	curved bristle	10mm	44.45mm	116.84mm	2 pcs
B	BR-8115	curved bristle	13.99mm	92.075mm	203.2mm	2 pcs
C	BR-8210	tapered fan bristle	4.75mm to 12.7mm	76.2mm	177.8mm	2 pcs
D	BR-8310	straight bristle	12.7mm	63.5mm	304.8mm	3 pcs
E	BR-8315	straight bristle	12.7mm	76.2mm	304.8mm	3 pcs
F	BR-8410	fan bristle	12.7mm	76.2mm	152.4mm	3 pcs


BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

BONE REAMER BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a rigid stainless steel handle with loop end.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-31-118	3mm	76mm	787mm	3 pcs
B	BR-31-197	5mm	76mm	787mm	3 pcs
C	BR-31-276	7mm	76mm	787mm	3 pcs




BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

POLYPROPYLENE BRUSHES

Channel cleaning brushes featuring rigid polypropylene bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	PP-16-039	1mm	25.4mm	406mm	2 pcs
	PP-22-039	1mm	25.4mm	558.8mm	2 pcs
B	PP-16-062	1.58mm	25.4mm	406mm	2 pcs
	PP-22-062	1.58mm	25.4mm	558.8mm	2 pcs
C	PP-16-087	2.2mm	25.4mm	406mm	2 pcs
	PP-22-087	2.2mm	25.4mm	558.8mm	2 pcs
D	PP-16-105	2.67mm	50.8mm	406mm	2 pcs
	PP-22-105	2.67mm	50.8mm	558.8mm	2 pcs
E	PP-16-125	3.18mm	50.8mm	406mm	2 pcs
	PP-22-125	3.18mm	50.8mm	558.8mm	2 pcs
F	PP-16-158	4mm	50.8mm	406mm	2 pcs
	PP-22-158	4mm	50.8mm	558.8mm	2 pcs
G	PP-18-197	5mm	50.8mm	457mm	2 pcs
	PP-22-197	5mm	50.8mm	558.8mm	2 pcs




BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

POLYPROPYLENE BRUSHES

Channel cleaning brushes featuring rigid polypropylene bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	PP-18-250	6.35mm	50.8mm	457mm	2 pcs
	PP-24-250	6.35mm	50.8mm	610mm	2 pcs
B	PP-18-393	10mm	50.8mm	457.2mm	2 pcs
C	PP-18-591	15mm	50.8mm	457.2mm	2 pcs
D	PP-18-750	19.05mm	152.4mm	457.2mm	2 pcs
E	PT-28-118 (pull through)	3mm	50.8mm	711mm	2 pcs
F	PT-28-197 (pull through)	5mm	50.8mm	711mm	2 pcs



A



B



C



D



E


(PULL THROUGH)



F

(PULL THROUGH)

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



Instrument Cleaning – Pipe Cleaners

PIPE CLEANERS - REEL


Available in a range of different materials. Rolls allow you to cut to the size you need, minimising waste. Single use.

	Product code	Hardness	Material	Diameter	Length	Qty
A	26322	soft	cotton	3mm	10000mm	1 reel
B	26323	medium	cotton/nylon	3mm	10000mm	1 reel
C	26328	medium	polyester/nylon	15mm	5000mm	1 reel
D	26329	medium	polyester/nylon	15mm	5000mm	1 reel
E	26334	hard	polyester/nylon	3mm	7500mm	1 reel
F	26335	hard	polyester/nylon	6mm	7500mm	1 reel
G	26336	hard	polyester/nylon	3mm	7500mm	1 reel
H	26337	medium	polyester/nylon	6mm	7500mm	1 reel
I	PC-300P	N/A	polyester	3.18mm	15240mm	1 reel
J	PC-400P	N/A	polyester	6.35mm	15240mm	1 reel

PIPE CLEANERS - LOOSE

Available in either polyester (non-woven) or polyester and nylon. Single use.

	Product code	Hardness	Material	Diameter	Length	Qty
K	26324	soft/hard	polyester/nylon	3mm	150mm	100 pcs
L	26325	soft	polyester (non-woven)	9mm	300mm	50 pcs
M	26326	soft	polyester (non-woven)	12mm	300mm	50 pcs
N	26327	medium	polyester/nylon	15mm	300mm	20 pcs
O	26342	soft	polyester (non-woven)	3mm	300mm	50 pcs



FLEXISTEM AND MICROCLEAN

Flexistem and Microclene are unique cleaners for all small-bore cannula. Various diameters from 3mm to 15mm are available constructed from very soft to hard bristles inserted between two soft mild steel twisted wires which assist in following contours while resisting buckling. Single use. Flexistem FLEX007A and FLEX007B has a wavy bristle configuration to reduce buckling which can occur in long instruments.

FLEXISTEM

Product code	Type	Material	Diameter	Length	Qty
A FLEX001	soft cleaner	cotton	3mm	10000mm	1 reel
B FLEX002	super bristle	polyester/nylon	3mm	10000mm	1 reel
FLEX003	double ended	polyester/nylon	3mm	300mm	100 pcs
FLEX004	fleece	polyester	9mm	300mm	50 pcs
C FLEX005	fleece	polyester	12mm	300mm	50 pcs
D FLEX006	bristle stem	polyester	15mm	300mm	20 pcs
E FLEX007	bristle stem	polyester/nylon	15mm	5000mm	1 reel
FLEX007A	bristle stem	polyester/nylon	15mm	5000mm	1 reel
FLEX007B	bristle stem	polyester/nylon	N/A	66.5mm	10 pcs
F FLEX033	fleece	polyester	3mm	300mm	50 pcs

MICROCLEAN


Product code	Type	Material	Diameter	Length	Qty
G FLEX021	super fine	polyester/nylon	3mm	7500mm	1 reel
H FLEX022	medium fine	polyester/nylon	6mm	7500mm	1 reel
I FLEX023	super brush	polyester/nylon	3mm	7500mm	1 reel
J FLEX024	medium brush	polyester/nylon	6mm	7500mm	1 reel

Instrument Cleaning – Toothbrush-Style Cleaning Brushes

TOOTHBRUSH-STYLE CLEANING BRUSHES

Antimicrobial nylon, brass, or stainless bristles on a sturdy, plastic handle. Designed for external cleaning.

Product code	Description	Total length	Qty
A	09320 extra strong synthetic bristles	165mm	10 pcs
B	09325 stainless steel bristles	180mm	1 pc
C	09485 extra strong synthetic bristles, angled handle	235mm	2 pcs
D	09478 soft synthetic bristles	215mm	1 pc
E	09488 synthetic bristles	180mm	1 pc
F	09489 brass bristles	180mm	1 pc
G	26491 soft synthetic bristles	215mm	1 pc
H	26492 hard synthetic bristles	215mm	1 pc
I	26029 nylon bristles, with sliding cover	55mm	3 pcs





A
B
C
D
E
F
G
H
I

Instrument Cleaning – Toothbrush-Style Cleaning Brushes

TOOTHBRUSH-STYLE CLEANING BRUSHES

Antimicrobial nylon, brass, or stainless bristles on a sturdy, plastic handle. Designed for external cleaning.

Product code	Description	Total length	Qty
A 09084	soft bristles, double-ended	175mm	2 pcs
B 09459	extra strong synthetic bristles, double-ended	175mm	2 pcs
C 09098	nylon bristles, double-ended	175mm	2 pcs
D 09460	extra strong, synthetic bristles	215mm	3 pcs
E B-2100	brass bristles, three rows	178mm	3 pcs
F N-2121	nylon bristles, double-ended, three rows/one row	178mm	2 / 50 pcs
G N-2221	nylon bristles, double-ended, three rows/one row	175mm	2 / 50pcs
H S-2121	stainless steel bristles, double-ended, three rows/one row	178mm	2 pcs
I S-2100	stainless steel bristles, three rows	178mm	3 pcs
J D-2100BN	brass, nylon bristles, three rows	178mm	3 pcs
K D-2100SN	stainless steel bristles, nylon, three rows	178mm	3 pcs



J
K



Instrument Cleaning – Toothbrush-Style and General Instrument Cleaning Brushes



A



B



C



D



E



F

TOOTHBRUSH-STYLE CLEANING BRUSHES

Antimicrobial nylon, brass, or stainless bristles on a sturdy, plastic handle. Designed for external cleaning.

Product code	Description	Total length	Qty
A N-2100	nylon bristles, three rows	178mm	3 pcs
B N-2200	extra rigid nylon bristles, angled handle	235mm	2 pcs
C S-2300	soft stainless steel bristles, three rows	178mm	3 / 50 pcs
D S-2131	soft stainless steel bristles, one row	178mm	3 / 50 pcs
79325	stainless steel bristles	180mm	1 pc
96281.3	stainless steel bristles	180mm	3 pcs
96281.4	extra strong synthetic bristles, angled handle	235mm	3 pcs

GENERAL INSTRUMENT CLEANING BRUSHES

A range of instrument cleaning brushes with strong durable nylon bristles and handles.

Product code	Description	Total length	Qty
E BBR001	nylon bristles	270mm	1 pc
F BBR002	nylon bristles	395mm	1 pc
G BBR003	nylon bristles	184mm	1 pc



G



Instrument Cleaning – General Instrument and Steriliser Cleaning Brushes

GENERAL INSTRUMENT CLEANING BRUSHES

Antimicrobial nylon or stainless steel bristles on a sturdy, plastic handle. Designed for external cleaning.

Product code	Description	Length	Qty
A	BR-2300 stainless steel	152.4mm	2 pcs
B	BR-04-003 stainless steel	91.44mm	2 pcs
C	09099 stainless steel with nylon bristles	450mm	2 pcs

STERILISER CLEANING BRUSH

Brush allows for greater reach and cleaning ability. Bristles are at a forward angle to loosen soil without marring surface.

Product code	Description	Bristle area	Length	Qty
D	BR-5000 double-surface, polyester bristles	152.4mm x 127mm	787.4mm	1 pc
E	BR-5200 threaded handle, plastic	N/A	914.4mm	1 pc

PISTOL-GRIP BRUSH

Soak-proof design with pistol-grip handle. Crimped nylon bristles on a block handle that resists chipping and cracking.

Product code	Description	Bristle length	Bristle area	Length	Qty
F	BR-4000 crimped nylon bristle, plastic pistol grip	50.8mm	101.6mm x 101.6mm	203.2mm	1 pc

STERILISER DRAIN CLEANING BRUSH

White nylon bristles designed to clean out drains, lines and openings.

Product code	Diameter	Diameter	Length	Qty
G	BR-5500 long handed brush, nylon bristles, with plastic handle	34.8mm	711.2mm	1 pc
H	BR-5600 drainline brush, nylon bristles, with stainless steel handle	50.8mm	609.6mm	1 pc


Instrument Cleaning – Brush Assortments

K

Instrument Cleaning




A




B



C



D



E



ORIGINAL BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- A BR-4199 BR-18-062, BR-18-087, BR-18-105, BR-18-125, BR-18-250, FT-18-500, N-2100, N-2121, S-2100

SMALL DIAMETER CHANNEL BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- B BR-4299 BR-08-158, BR-12-079, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-24-197, BR-24-200

LARGE DIAMETER CHANNEL BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- C BR-4399 BR-06-375, BR-08-630, BR-16-276, BR-16-393, BR-18-250, BR-18-591, BR-18-880, BR-24-312, BR-24-440

POPULAR CHANNEL BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- D BR-4499 BR-06-375, BR-08-158, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-18-250, BR-24-200

MICRO-OPHTHALMIC BRUSH ASSORTMENT






Product code Includes one of each of the following brushes

- E BR-4599 BR-18-039, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-18-158, N-2100

STERILISER CLEANING BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- BR-6099 BR-5000, BR-5200, BR-5400, BR-5500, BR-5600



Instrument Cleaning – Endoscope Cleaning Brushes

ENDOSCOPE BRUSHES

Nylon bristles attached to flexible, coloured plastic tubing. Single-ended. Single use. Available in 1 pc or 50 pcs.

Product code	Colour	Brush diameter	Bristle area length	Length	Qty	
A	BR-075-18	blue	1.8mm	20mm	750mm	1 pc
B	BR-075-24	yellow	2.4mm	20mm	750mm	1 pc
C	BR-075-32	green	3.2mm	20mm	750mm	1 pc
D	BR-075-42	pink	4.2mm	20mm	750mm	1 pc
E	BR-100-18	blue	1.8mm	20mm	1000mm	1 pc
F	BR-100-24	yellow	2.4mm	20mm	1000mm	1 pc
G	BR-100-32	green	3.2mm	20mm	1000mm	1 pc
H	BR-100-42	pink	4.2mm	20mm	1000mm	1 pc
I	BR-150-18	blue	1.8mm	20mm	1500mm	1 pc
J	BR-150-24	yellow	2.4mm	20mm	1500mm	1 pc
K	BR-150-32	green	3.2mm	20mm	1500mm	1 pc
L	BR-150-42	pink	4.2mm	20mm	1500mm	1 pc
M	BR-200-18	blue	1.8mm	20mm	2000mm	1 pc
N	BR-200-24	yellow	2.4mm	20mm	2000mm	1 pc
O	BR-200-32	green	3.2mm	20mm	2000mm	1 pc
P	BR-200-42	pink	4.2mm	20mm	2000mm	1 pc

ENDOSCOPE BRUSHES

Tapered cleaning brushes. Nylon bristles attached to flexible plastic tubing. Single-ended. Single use. Available in 1 pc and 50 pcs.








Product code	Brush diameter	Bristle area length	Length	Qty	
Q	BR-3308	tapered 3mm to 5mm	10mm	2400mm	1 pc
R	BR-3309	tapered 3mm to 5mm	10mm	270mm	1 pc
S	BR-3310	tapered 3mm to 5mm	20mm	2400mm	1 pc

Instrument Cleaning – Endoscope Cleaning Brushes

ENDOSCOPE BRUSHES

A variety of single, double, and triple bristled endoscope brushes as well as valve brushes and kits.

Product code	Description	Brush diameter	Bristle area length	Length	Qty
A	EB10100	double-ended	10mm & 5mm	30mm & 20mm	2500mm 50 pcs
B	EB10102	double-ended	3mm	20mm	2300mm 50 pcs
C	EB10103	double-ended	6mm	20mm	2500mm 50 pcs
D	EB10103.3	triple-cleaning	6mm	20mm	2500mm 50 pcs
E	EB10104	double-ended with handle	10mm & 5mm	40mm & 20mm	2300mm 25 pcs
F	EB10105	double-ended with pipe cleaner	5mm	125mm & 20mm	2300mm 50 pcs
G	EB10106	channel brush with handle	2mm	10mm	1200mm 50 pcs
H	EB10136	double-ended	3mm	20mm	1200mm 50 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Endoscope Cleaning Brushes

ENDOSCOPE BRUSHES




A variety of endoscope brushes conveniently packaged with a valve brush. Single use.

Product code	Description	Qty
A EB10103.3C	includes triple-cleaning channel EB10103.3 and valve brush EB10101	50 pcs
B EB10107C	includes double-ended channel EB10103 and valve brush EB10101	50 pcs
C EB10107G	includes double-ended channel EB10102 and valve brush EB10101	50 pcs


ENDOSCOPE BRUSHES (TAPERED, SINGLE-ENDED)

Nylon bristles attached to flexible plastic tubing. Single-ended. Single use.


Product code	Description	Bristle area length	Length	Qty
BR-3711-50	tapered 4mm to 6mm	10.5mm	2400mm	50 pcs
D BR-3799-50	BR-3711-50 and BR-3201	N/A	2400mm	50 pcs



BRUSH TIP DETAILS
SHOWN ACTUAL SIZE



Instrument Cleaning – Endoscope Cleaning Brushes



ENDOSCOPE BRUSHES (TAPERED, DOUBLE-ENDED)

Nylon bristles attached to flexible plastic tubing. Double-ended. Single use. Available in 1 pc or 50 pcs.

Product code	Brush diameter	Bristle area length	Length	Qty
A BR-3320	tapered 3mm to 5mm	10mm	2400mm	1 pc
B BR-3330	tapered 3mm to 5mm	20mm	2400mm	1 pc

ENDOSCOPE BRUSHES (NON-TAPERED, DOUBLE-ENDED)


Nylon bristles attached to flexible plastic tubing. 3mm brush on one end, 6mm brush on the other. Single use. Available in 1 pc or 50 pcs.

Product code	Brush diameter	Bristle area length	Length	Qty
C BR-3340	3mm & 6mm	20mm	2400mm	1 pc

PULL-THROUGH ENDOSCOPE BRUSHES

Delicate pull-through cleaning brushes. Nylon bristles attached to twisted stainless steel wire. Single-ended. Single use. Available in 1 pc or 50 pcs.

Product code	Brush diameter	Bristle area length	Length	Qty
D BR-3030	1.1mm	12.7mm	650mm	1 pc
BR-3035	1.1mm	10.0mm	2300mm	1 pc
E BR-3040	2.0mm	12.7mm	1200mm	1 pc
F BR-3410	3.0mm	13.0mm	1200mm	1 pc
G BR-3420	3.0mm	13.0mm	1600mm	1 pc



Instrument Cleaning – Endoscope Cleaning Brushes



ALBARRAN BRUSH

Delicate brush designed for cleaning elevator channel on a duodenoscope. Single use.

Product code	Description	Brush diameter	Bristle area length	Qty
A EB10108	albarran entrance brush	2.3mm	13mm	100 pcs

VALVE BRUSH





Short brush designed for cleaning valve ports and suction wells on an endoscope. Single use.

Product code	Description	Brush diameter	Bristle area length	Qty
B BR-3201	suction well brush	11mm	45mm	3 / 50 pcs
C BR-3202	double-ended	4.3mm & 11mm	10.2mm & 45mm	3 / 50 pcs
D EB10101	double-ended	10mm & 5mm	30mm & 20mm	100 pcs
E BR-3307	finger loop style	11mm	45mm	1 pc

DID YOU KNOW?

Key Surgical is a global leading partner to hospitals and theatres or procedural facilities, supporting the needs in endoscopy, operating room, and instrument reprocessing. Our comprehensive portfolio is designed with a focus on patient safety; allowing our Customers to efficiently source the high quality products they need with the best Customer service and overall experience with every interaction. Our UK office is based in High Wycombe but our main headquarters are in the US and our products are distributed in over 70 countries around the world. Throughout our continued growth we have remained true to our mission; to be the easiest company in our industry with which to do business.

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



Instrument Cleaning – Endoscope Cleaning Brushes

PREMIUM ENDOSCOPE BRUSH

Nylon bristles on twisted stainless steel wire or spring coil for hard-to-turn corners. Loop handle at one end and stainless steel beaded tip at bristle end. Reusable.

Product code	Description	Brush diameter	Bristle area length	Qty
A BR-3300	spring coil	tapered 4mm to 6mm	10mm	1 pc
B BR-3305	spring coil	tapered 3mm to 5mm	10mm	1 pc
C BR-3702	spring coil	6mm	5mm	1 pc
D BR-3701	twisted wire	2.4mm	6mm	1 pc
E BR-3703	twisted wire	2.4mm	6mm	1 pc

ENDOSCOPIC BRUSH ASSORTMENTS

Product code	Includes one of each of the following brushes
BR-4799	BR-12-079, BR-12-118, BR-12-197, BR-16-197, BR-16-393, BR-24-197, BR-24-393, BR-3202, N-2121
F BR-100-99	BR-100-18, BR-100-24, BR-100-32, BR-100-42
G BR-150-99	BR-150-18, BR-150-24, BR-150-32, BR-150-42
H BR-200-99	BR-200-18, BR-200-24, BR-200-32, BR-200-42

