

STERILISATION TRAYS AND ACCESSORIES

Our range of high-quality stainless steel washing and instrument trays and baskets are produced in line with industry standard sizes and can be made to fit customer requirements subject to a minimum order quantity.

The sterilisation tray range holds smaller items that are difficult to keep track of during sterilisation as well as larger items with small parts. Our mesh trays and instrument trays are heavy duty stainless steel for tough, clean, and corrosion-resistant surfaces.

Nic-Nac Baskets are designed to transport small items through the wash and autoclave process. All Nic-Nac Baskets are produced with a perforated stainless-steel sheet with a scalloped hinged lid.

CBM basket ranges are manufactured with AIS1304 grade stainless steel and can be used within sterilisation containers or with sterilisation wrap. Due to electro polishing, they are resistant to colour changes and have a surface treatment that increases their water drainage capacity.

Plastic Sterilisation Trays are designed to cover most requirements for protecting instruments during sterilisation and transportation. Grid holes in the tray provide drainage, fixing positions for various inserts, and buttons to raise silicone pin mats from the base. Transparent lids are strongly fitted to provide additional security and ensure all contents are highly visible. Suitable for low temperature sterilisation processes.

We also have our range of accessories such as Corner Protectors, Protective Mats, Soaker Sheets and Labels that are effective solution to reduce damage of trays but also instruments.

Trinket Baskets

Nic-Nac Baskets

Sieve Trays

Flat Base Baskets

Mesh Baskets

Micro Mesh Trays

Stainless Steel Baskets

CBM Baskets

Plastic Sterilisation Trays

Stainless Steel Scope Baskets

Silicone Inserts

Corner Protectors

Soaker Sheets

Protective Mats

Labels

Please note: Unless specified, all products are not made with natural rubber latex.



Sterilisation Trays and Accessories - Trinket Baskets



TRINKET BASKETS

Small, stainless steel trinket baskets to hold fine items such as screws and bolts. Available with lids and designed to transport small items through the washing and autoclave process. Compatible with hydrogen peroxide, EtO and steam sterilisation up to 137°C.

Product code	Description	Size	Qty
A CIL0705	trinket box	90 x 60 x 25mm	1 pc
CIL0706	trinket box, perforated lid	105 x 70 x 25mm	1 pc
CIL0707	trinket box	120 x 80 x 25mm	1 pc
B CIL76-0175	trinket box, perforated sides	70 x 30 x 15mm	1 pc
CIL76-0185	trinket box, wire mesh lid	90 x 60 x 25mm	1 pc
CIL76-0186	trinket box, wire mesh lid	105 x 70 x 25mm	1 pc
CIL76-0187	trinket box, wire mesh lid	120 x 80 x 25mm	1 pc
CIL76-0189	trinket box, perforated sides and base	150 x 100 x 25mm	1 pc
C CIL0406	fine mesh basket	80 x 80 x 35mm	1 pc
D CIL76-0204	fine mesh basket	105 x 70 x 25mm	1 pc
SJ/165-009	fine mesh basket, with lid	80 x 80 x 60mm	1 pc



Sterilisation Trays and Accessories - Nic-Nac Baskets

NIC-NAC BASKETS

A range of small Nic-Nac Baskets manufactured from fine perforated stainless steel with scalloped hinged lids and catch. Designed to accommodate those small items that get lost in a tray and for transportation through the washing process. Compatible with hydrogen peroxide, EtO and steam sterilisation up to 134°C.

Product code	Size	Qty
A NIC001R	80 x 30mm	1 pc
B NIC002	178 x 127 x 25mm	1 pc
C NIC003	229 x 178 x 25mm	1 pc
C NIC004DEEP	180 x 180 x 50mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays




SIEVE FOR SMALL PARTS

Small part sieves and accessories to help you to organise your work processes. Can be used for reprocessing, transportation, and storage for the protection of instruments or small parts. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
A 08405	without lid and clips	236 x 121 x 46mm	1 pc
08947	lid for ref 08405, with clips	236 x 121 x 8mm	1 pc
B 08406	set with lid and clips	121 x 121 x 29mm	1 pc
C 08408	set with lid and clips	121 x 73 x 29mm	1 pc
D 08409	without lid and clips	121 x 73 x 29mm	1 pc
E 08430	lid for ref 08409, with clips	121 x 73 x 7mm	1 pc
08448	set with lid and clips	121 x 121 x 29mm	1 pc
F 08449	set with lid and clips	121 x 121 x 47mm	1 pc
G 08450	without lid and clips	121 x 121 x 29mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays



SIEVE FOR SMALL PARTS

Small part sieves and accessories to help you to organise your work processes. Can use for reprocessing, transportation, and storage for the protection of instruments or small parts. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
08452	lid for ref 08450, with clips	121 x 121 x 7mm	1 pc
08451	without lid and clips	121 x 121 x 47mm	1 pc
08454	with lid and clips	236 x 73.5 x 28mm	1 pc
A 08457	without lid and clips	236 x 73,5 x 28mm	1 pc
08461	lid for ref 08457, with 2 clips	236 x 73.5 x 7mm	1 pc
B 08463	set with lid and clips	236 x 121 x 28mm	1 pc
C 08464	set with lid and clips	236 x 121 x 28mm	1 pc
08465	without lid and clips	236 x 121 x 28mm	1 pc
08467	lid ref 08465, with 2 clips	236 x 121 x 7mm	1 pc
08466	without lid and clips	236 x 121 x 28mm	1 pc
08468	lid for ref 08466, with 2 clips	236 x 121 x 7mm	1 pc
D 08442	set with lid and clips	63 x 63 x 23.4mm	1 pc
E 08434	with lid	80 x 80 x 40mm	1 pc
F 08435	with lid	80 x 40 x 30mm	1 pc
G 08896	with lid	40 x 40 x 25mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays

SIEVE FOR SMALL PARTS

Small part sieves and accessories to help you to organise your work processes. Can be used for reprocessing, transportation, and storage for the protection of instruments or small parts. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
A 08898	with lid	160 x 40 x 20mm	1 pc
B 08899	with lid	200 x 40 x 40mm	1 pc
C 08900	with lid	80 x 40 x 40mm	1 pc
D 08901	with lid	80 x 80 x 40mm	1 pc
E 08906	with lid	250 x 200 x 25mm	1 pc
08907	with lid	300 x 180 x 25mm	1 pc
F 08772	with lid	210 x 80 x 40mm	1 pc
08948	with lid	200 x 60 x 40mm	1 pc
08875	with lid	300 x 150 x 70mm	1 pc
G 08950	with lid	300 x 150 x 70mm	1 pc
H 08762	without lid	200 x 50 x 50mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays



SIEVE FOR SMALL PARTS WITH LENGTH DIVIDERS


Small part sieves and accessories to help you to organise your work processes. Can be used for reprocessing, transportation, and storage for the protection of instruments or small parts. Lengthwise dividers for small parts sieves. For individual division. Includes two clips.

Product code	Description	Size	Qty
A 08410	lengthwise divider 1/1	224 x 24mm	1 pc
08412	lengthwise divider 1/1	224 x 42mm	1 pc
B 08482	lengthwise divider 2/3	300 x 43mm	1 pc
C 08436	lengthwise divider 1/1	465 x 43mm	1 pc
D 08414	crosswise divider 1/1	109 x 24mm	1 pc
E 08437	crosswise divider 1/3	71 x 43mm	1 pc
F 08481	crosswise divider 1/1	61.5 x 24mm	1 pc
08483	crosswise divider 1/1	232 x 43mm	1 pc
08484	crosswise divider 2/3	152 x 43mm	1 pc
08485	crosswise divider 2/3	152 x 43mm	1 pc

COMPARTMENTAL SIEVE FOR SMALL PARTS

A Compartmental Sieve Tray for small parts with 65 inner compartments. Made from stainless steel.

Product code	Description	Size	Qty
G 08880	65 inner compartments, with lid, hinge, spring cap and side walls	40 x 15mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays



SIEVE TRAYS

Our range of stainless steel Sieve Trays, with and without lids, helps you guarantee optimal and efficiency in your work process. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
A 08808	without lid	480 x 250 x 100mm	1 pc
08809	without lid	560 x 250 x 50mm	1 pc
08810	without lid	560 x 250 x 70mm	1 pc
B 08415	without lid	240 x 250 x 50mm	1 pc
08416	without lid	480 x 250 x 50mm	1 pc
08417	without lid	480 x 250 x 70mm	1 pc
08580	with lid	540 x 250 x 50mm	1 pc
C 08761	with lid	400 x 280 x 70mm	1 pc
08581	with lid	540 x 250 x 100mm	1 pc
D 08801	rochard sieve, with lid	725 x 275 x 113mm	1 pc





Sterilisation Trays and Accessories - Sieve Trays

SIEVE TRAYS AND MATCHING LIDS

Sieve Trays with matching lids, made of Chrome Nickel Steel including accessories as well. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
A 08455	set 1/6 XL, with lid, set comprising: 1 x system small parts sieve 1/6 XL with lid and clips	396 x 73.5 x 47mm	1 pc
B 08456	1/6, without lid	396 x 73.5 x 47mm	1 pc
C 08459	1/6 XL, without lid	396 x 73.5 x 47mm	1 pc
D 08453	set 1/6, with lid, set comprising: 1 x system small parts sieve 1/6 with lid and clips	396 x 73.5 x 47mm	1 pc
E 08460	lid for ref 08456 with 2 clips	236 x 73.5 x 7mm	1 pc
08462	lid for ref 08459 with 2 clips	396 x 73.5 x 7mm	1 pc
08470	set 1/3L, with lid, set comprising: 1 x system small parts sieve 1/3L with lid and clips	466 x 73.5 x 47mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays

SIEVE TRAY SETS

Sieve Trays Sets, made of Chrome Nickel Steel including accessories as well. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
A 08469	1 x system small parts sieve with lid and clips	466 x 73.5 x 28mm	1 pc
B 08471	1 x system small parts sieve with lid and clips	236 x 236 x 28mm	1 pc
C 08472	1 x system small parts sieve with lid and clips	236 x 236 x 46mm	1 pc

SIEVE TRAYS WITH DROP HANDLES

Sieve Trays with drop handles are made of stainless steel and without lids. Compatible with steam sterilisation up to 134°C.

Product code	Size	Qty
D 08908	400 x 250 x 70mm	1 pc
E 08909	250 x 240 x 50mm	1 pc



Sterilisation Trays and Accessories - Flat Base Baskets



FLAT BASE BASKETS

Our Flat Base Baskets range are stainless steel for instrument washing, disinfection, sterilisation, storage, and handling. Compatible with steam sterilisation processes up to 137°C.

Product code	Description	Size	Qty
CIL76-0206	flat base basket, perforated sides	150 x 50 x 35mm	1 pc
A CIL76-0049	flat base basket, perforated sides	480 x 250 x 70mm	1 pc
CIL76-0039	flat base basket, perforated sides	380 x 240 x 50mm	1 pc
CIL76-0047	flat base basket, perforated sides	480 x 250 x 50mm	1 pc
CIL76-0048	flat base basket, perforated sides	480 x 250 x 60mm	1 pc
B CIL76-0030	flat base basket, perforated sides	250 x 240 x 50mm	1 pc
CIL76-0057	flat base basket, perforated sides	480 x 360 x 60mm	1 pc
CIL76-0058	flat base basket, perforated sides	480 x 480 x 50mm	1 pc
CIL76-0060	flat base basket, perforated sides	480 x 480 x 70mm	1 pc
CIL76-0066	flat base basket, perforated sides	520 x 250 x 50mm	1 pc
CIL76-0076	flat base basket, perforated sides	550 x 250 x 60mm	1 pc
CIL76-0079	flat base basket, perforated sides	550 x 430 x 50mm	1 pc
CIL76-0086	flat base basket, perforated sides	600 x 300 x 60mm	1 pc
CIL76-0087	flat base basket, perforated sides	600 x 350 x 60mm	1 pc
CIL76-0088	flat base basket, perforated sides	670 x 300 x 70mm	1 pc
CIL76-0089	flat base basket, perforated sides	700 x 250 x 60mm	1 pc
CIL76-0090	flat base basket, perforated sides	700 x 300 x 90mm	1 pc
CIL76-0095	lid for flat base basket	250 x 240mm	1 pc
CIL76-0096	lid for flat base basket	480 x 250mm	1 pc
CIL76-0097	lid for flat base basket	480 x 480mm	1 pc
CIL0202LID	lid for wire mesh basket	240 x 250mm	1 pc
CIL0210LID	lid for wire mesh basket	480 x 250mm	1 pc
C CIL0202	wire mesh basket	250 x 240 x 60mm	1 pc
D CIL0210	wire mesh basket	480 x 250 x 60mm	1 pc

FLAT BASE DIN BASKETS

Product code	Description	Size	Qty
SJ-131-020	din basket, with fitted lid on top	650 x 250 x 70mm	1 pc
SJ-275-275	stainless steel din basket perforated sides, flat base	275 x 275 x 50mm	1 pc



Sterilisation Trays and Accessories - Mesh Baskets



MESH BASKETS AND LIDS

Stainless steel Fine Mesh Baskets with lids for instrument washing, disinfection, sterilisation, storage, and handling. Compatible with steam sterilisation up to 137°C, and plasma sterilisation processes.

Product code	Description	Size	Qty
CIL76-0143	all wire mesh basket	480 x 250 x 50mm	1 pc
CIL76-0200	all wire mesh basket	80 x 40 x 30mm	1 pc
A CIL76-0202	fine mesh basket	80 x 80 x 55mm	1 pc
CIL76-0210	all wire mesh basket	190 x 100 x 85mm	1 pc
CIL76-0215	fine mesh basket	200 x 100 x 35mm	1 pc
CIL76-0222	fine mesh basket	240 x 90 x 40mm	1 pc
CIL76-0250	fine mesh basket	120 x 155 x 40mm	1 pc
CIL76-0253	fine mesh lid	120 x 155mm	1 pc
CIL76-0256	fine mesh basket, take off lid, hinges and locking clip	120 x 155 x 50mm	1 pc
CIL76-0252	fine mesh basket, with lid	235 x 155 x 40mm	1 pc
CIL76-0475	all wire mesh basket	250 x 170 x 70mm	1 pc


PERFORATED PLATE MESH TRAYS

Stainless steel Perforated Plate Mesh Trays. Suited for soft packaging. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
B 08432	with drop handles, without lid	480 x 250 x 105mm	1 pc
08473	with drop handles, without lid	480 x 250 x 60mm	1 pc
08474	with drop handles, without lid	480 x 250 x 65mm	1 pc
08475	with drop handles, without lid	480 x 250 x 100mm	1 pc
C 08476	with drop handles, without lid	240 x 250 x 60mm	1 pc
08478	with drop handles, without lid	240 x 250 x 65mm	1 pc
08579	with turn-lock and lid	478 x 247mm	1 pc
08583	without lid	480 x 250 x 100mm	1 pc
D 08934	flushing mesh basket, with lid	240 x 250 x 60mm	1 pc
E 08942	flushing mesh basket, with lid	480 x 250 x 60mm	1 pc
08943	flushing mesh basket, with lid	240 x 250 x 60mm	1 pc
08945	flushing mesh basket, with lid	220 x 140 x 54mm	1 pc
08975	without lid	120 x 250 x 60mm	1 pc
08976	without lid	240 x 250 x 66mm	1 pc



Sterilisation Trays and Accessories - Micro Mesh Trays



MICRO MESH TRAYS

Micro Mesh Trays are designed to hold smaller items during washing, sterilisation, storage, and transportation. Micro mesh allows you to keep larger items that have small loose parts together. Micro mesh size is approximately 2mm x 2mm. Compatible with steam and gas plasma sterilisation.

Product code	Description	Size	Qty
A MT-5500	with lids	40 x 40 x 20mm	1 pc
B MT-5600	with lids	80 x 40 x 20mm	1 pc
C MT-5650	with lids	160 x 40 x 20mm	1 pc
D MT-5700	with lids	70 x 70 x 30mm	1 pc
E MT-5800	with lids	80 x 80 x 40mm	1 pc
F MT-5900	with lids	90 x 90 x 50mm	1 pc
G MT-6020	with lids	180 x 130 x 20mm	1 pc
H MT-7100	with lids	225 x 110 x 50mm	1 pc
I MT-7200	with lids	270 x 135 x 40mm	1 pc
J MT-7300	with lids	400 x 250 x 50mm	1 pc
K MT-9001	with drop handles	480 x 250 x 50mm	1 pc
L MT-9101	with drop handles	250 x 240 x 50mm	1 pc
M MT-9521	with drop handles	700 x 210 x 50mm	1 pc



Sterilisation Trays and Accessories - Stainless Steel Baskets

PREMIUM STAINLESS STEEL BASKETS


Stainless Steel Baskets for holding instruments during sterilisation. Compatible with steam sterilisation up to 137°C.



Product code	Description	Size	Qty
A BSK316	stainless steel basket	480 x 250 x 50mm	1 pc
B BSK317	stainless steel basket	250 x 250 x 50mm	1 pc
C BSK318	stainless steel basket	480 x 480 x 50mm	1 pc
LID316	stainless steel perforated lid to fit BSK316	480 x 250mm	1 pc
LID317	stainless steel perforated lid to fit BSK317	250 x 250mm	1 pc
LID318	stainless steel perforated lid to fit BSK318	480 x 480mm	1 pc
BSK316H	stainless steel basket, with handle	480 x 250 x 50mm	1 pc
BSK317H	stainless steel basket, with handle	250 x 250 x 50mm	1 pc
BSK318H	stainless steel basket, with handle	480 x 480 x 50mm	1 pc
BELBSK116	stainless steel basket, mesh base with handles	480 x 250 x 60mm	1 pc
BELBSK117	stainless steel basket, with handle	250 x 240 x 60mm	1 pc



Sterilisation Trays and Accessories - Stainless Steel Baskets



STRINDBERG TRAYS

A series of stainless steel trays designed to protect delicate instruments during washing and sterilisation. Instruments are held in place with a firm clip protected with a silicone shield.

Product code	Description	Size	Qty
A RON005	full tray, 20 instrument slots	275 x 180 x 25mm	1 pc
RON006	half tray, 9 instrument slots	131 x 180 x 50mm	1 pc
RON007	quarter tray, 5 instrument slots	73 x 180 x 25mm	1 pc


TRAYS AND CASSETTES

Stainless steel tray for the storage of instruments.

Product code	Description	Size	Qty
B 08940	storage tray for instruments, with silicone holder 5-fold with fixing bracket	180 x 85 x 23mm	1 pc




Sterilisation Trays and Accessories - CBM Baskets



MINI AND MICRO BASKETS

Small, stainless steel baskets to hold small items such as screws and bolts. Available with lids and designed to transport small items through the washing and autoclave process.

	Product code	Description	Size	Qty
A	CBM63J	mini basket with lid, without handles	200 x 80 x 35mm	1 pc
B	CBM67F2	micro basket with lid, without handles	170 x 135 x 35mm	1 pc
C	CBM63M	mini basket with lid, without handles	220 x 60 x 35mm	1 pc
D	CBM76D	mini basket with lid and handles	267 x 125 x 49mm	1 pc



STERILISATION BASKETS SYSTEM

CBM Sterilisation Baskets Systems are manufactured with stainless steel and can be used within sterilisation containers or with sterilisation wrap. Due to electro polishing, they are resistant to colour changes and have a surface treatment that increases their water drainage capacity.

	Product code	Description	Size	Qty
E	CBM71/1	stainless steel basket, without lid	240 x 240 x 50mm	1 pc
	CBM71LID	stainless steel basket, with lid	240 x 240mm	1 pc
	CBM74/1	stainless steel basket, without lid	480 x 240 x 50mm	1 pc
	CBM74LID	stainless steel basket, with lid	480 x 240mm	1 pc
	CBM75/1	stainless steel basket, without lid	540 x 240 x 50mm	1 pc
	CBM75LID	stainless steel basket, with lid	540 x 240mm	1 pc
	CBM79/2	stainless steel basket, without lid	650 x 240 x 70mm	1 pc
	CBM79/LID	stainless steel basket, with lid	650 x 240mm	1 pc

SCOPE BASKETS SYSTEM

A range of baskets for the reprocessing, storage, and transportation of arthroscopes and endoscopes, available with integrated silicone instrument holders.

	Product code	Description	Size	Qty
F	CBM48/S	scope basket for 1 scope, with lid	480 x 150 x 80mm	1 pc
G	CBM66/S	scope basket for 1 scope, with lid	450 x 70 x 60mm	1 pc
H	CBM64ECV	empty scope basket with lid and handles	650 x 150 x 70mm	1 pc
I	CBM66AV	empty scope basket with lid and handles	580 x 150 x 80mm	1 pc
	CBM66/FEN	flexible scope basket with lid and 4 silicone supports	580 x 150 x 80mm	1 pc



Sterilisation Trays and Accessories - Plastic Sterilisation Trays

STERILISATION TRAYS

Sterilisation trays made of plastic and stainless steel. The plastic trays and cassettes are used for sterilisation and instrument storage.

Product code	Description	Colour	Size	Qty
A 08422	lid black-transparent with 4mm holes	beige	150 x 65 x 26mm	1 pc
B 08424	lid black-transparent with 6mm holes	beige	230 x 175 x 42mm	1 pc
C 08425	lid black-transparent with 6mm holes	beige	300 x 200 x 50mm	1 pc

WASH TRAYS

Wash Trays made for Da Vinci® systems.

Product code	Description	Size	Qty
VAS4XI	washing tray for Endo Wrist®	660 x 212 x 50mm	1 pc
08798	sieve tray	640 x 340 x 100mm	1 pc
08132	mesh basket with lid, made of stainless steel	670 x 100 x 70mm	1 pc



Sterilisation Trays and Accessories - Stainless Steel Scope Baskets


STAINLESS STEEL SCOPE BASKETS

Stainless steel wire mesh baskets for one or two endoscopes or arthroscopes with integral silicone instrument holders. This lockable range protects scope instruments whilst washing and sterilising.

Product code	Description	Size	Qty
A SB-1210	with silicone inserts	180 x 80 x 52mm	1 pc
B SB-4210	with silicone inserts	460 x 80 x 52mm	1 pc
C SB-4610	with silicone inserts	460 x 160 x 52mm	1 pc
D SB-6520	with silicone inserts, includes small lens basket with fixed bracket	640 x 150 x 77mm (80 x 80 x 40mm)	1 pc
CIL2728	scope basket	290 x 80 x 52mm	1 pc
CIL2727	scope basket	460 x 80 x 52mm	1 pc
CIL2726	scope basket	670 x 80 x 52mm	1 pc
CIL2701	scope basket, with small mesh basket	640 x 150 x 77mm	1 pc
CIL2729	scope basket for 2 endoscopes	460 x 80 x 52mm	1 pc
CIL2730	double scope basket	290 x 80 x 52mm	1 pc
CIL2751	double scope basket	290 x 160 x 52mm	1 pc
CIL2753	double scope basket	460 x 160 x 52mm	1 pc
CIL76-0001	endoscope / arthroscope basket	640 x 150 x 77mm	1 pc
CIL76-0010	universal scope basket	580 x 160 x 52mm	1 pc
CIL76-0011	universal scope basket	290 x 160 x 52mm	1 pc
CIL76-0012	universal scope basket	460 x 160 x 52mm	1 pc
SJ-102-001P	stainless steel, with 1 pin	460 x 80 x 52mm	1 pc
SJ-102-002P	stainless steel, with 2 pin	460 x 80 x 52mm	1 pc
SJ-103-002	with small mesh basket	640 x 150 x 77mm	1 pc



Sterilisation Trays and Accessories - Silicone Inserts



SILICONE INSERTS


For custom configuration within a tray to accommodate various sizes of instruments.

Product code	Description	Opening diameter	Size	Qty
A A237593	straight opening, 3 slots	N/A	57 x 15mm	1 pc
B A237596	straight opening, 6 slots	N/A	57 x 15mm	1 pc
A387596	straight opening, 6 slots	N/A	95 x 19mm	1 pc
A3875912	straight opening, 12 slots	N/A	95 x 19mm	1 pc
C A237591K	round opening, 1 slot	6mm	57 x 15mm	1 pc
D A237592K	round opening, 2 slots	6mm	57 x 15mm	1 pc
E A387592	round opening, 2 slots	6mm	95 x 19mm	1 pc
A2616902	round opening, 2 slots	10mm	63 x 25mm	1 pc
A2616901	round opening, 2 slots	10 & 13mm	53 x 25mm	1 pc
A387593	round opening, 3 slots	6mm	95 x 19mm	1 pc
F A237593K	round opening, 3 slots	6mm	57 x 15mm	1 pc
A387594	round opening, 4 slots	6mm	95 x 19mm	1 pc
A387595	round opening, 5 slots	6mm	95 x 19mm	1 pc
G A32759P2	long, round opening, 1 slot	7mm	89 x 32mm	1 pc
H A32759P1	long, round opening, 1 slot	16mm	89 x 32mm	1 pc
I A327593	long, round opening, 3 slots	7mm	89 x 32mm	1 pc
J A327593L	long, round opening, 3 slots	16 & 7 & 13mm	89 x 32mm	1 pc



Sterilisation Trays and Accessories - Corner Protectors

A



SURGICAL TRAY CORNER PROTECTORS

Tray Corner Protectors are made from medical-grade silicone and are intended for use in the organization and storage of surgical trays. Compatible with steam sterilisation. Single use.

Product code	Description	Colour	Size	Qty
A CS-9125C	foam, 4 pre-cut corners	blue	152 x 152 x 32mm	80 sheets

CORNER PROTECTORS

Tray Corner Protectors are made from medical-grade silicone and are intended for use in the organization and storage of surgical trays. Heat-resistant up to 134°C.

Product code	Material	Colour	Size	Qty
B 42006	silicone	pink	73.1 x 73.1 x 18.5mm	80 pcs
42007	card	white	100 x 100 x 50mm	600 pcs
C 42008	foam	blue	75 x 75 x 32mm	288 pcs
42010	card	white	100 x 100 x 50mm	600 pcs
42011	silicone	yellow	76 x 54 x 33mm	4 pcs
D 50062	silicone	violet	76 x 54 x 33mm	4 pcs

CARD AND CREPE CORNER PROTECTORS

Tray Corner Protectors are intended for use in the organization and storage of surgical trays. Suitable for steam sterilisation, they are available in crepe or card.

Product code	Material	GSM	Qty
CC001	card	171gsm	200 pcs
E CC002	crepe	171gsm	100 pcs

B




C



D



E



A



Sterilisation Trays and Accessories - Soaker Sheets

SOAKER SHEETS



140gsm Soaker Sheets are for use under trays to protect inner wraps from sharp corners and feet on trays and to absorb condensate during sterilisation. Soaker Sheets can also be used as a shelf liner for added protection. Single use.



Product code	Size	Qty
A SOK001SS	750 x 500mm	50 sheets
SOK002SS	500 x 375mm	100 sheets
SOK003SS	610 x 305mm	100 sheets
SOK004SS	500 x 200mm	100 sheets
SOK005SS	500 x 500mm	100 sheets
SOK006SS	250 x 250mm	100 sheets
SOK007SS	500 x 280mm	100 sheets
SOK008SS	1219 x 560mm	50 sheets
SOK009SS	1524 x 560mm	50 sheets
SOK010SS	914 x 508mm	50 sheets
SOK011SS	914 x 588mm	50 sheets
SOK012SS	1000 x 600mm	50 sheets
SOK013SS	1295 x 350mm	50 sheets
SOK016SS	200 x 125mm	1000 sheets



Sterilisation Trays and Accessories - Protective Mats



ABSORBENT PROTECTIVE MATS

Absorbent Mats are intended for placement on the bottom of a sterilisation tray to absorb wetness and humidity during processing. The 'dry' version, with an added absorbent layer, is also available to provide rapid moisture dispersion during steam sterilisation. Compatible with steam sterilisation. Single use.

PURPLE

Product code	Description	Size	Qty
A MAT001	double din	750 x 500mm	50 sheets
MAT002	double din	500 x 375mm	100 sheets
MAT003	double din	610 x 305mm	100 sheets
MAT004	din	500 x 250mm	100 sheets
MAT005	half din	250 x 250mm	100 sheets
MAT007	double din	500 x 500mm	100 sheets




PURPLE DRY

Product code	Description	Size	Qty
A DRYMAT001	double din	750 x 500mm	50 sheets
DRYMAT002	double din	500 x 375mm	100 sheets
DRYMAT003	double din	610 x 305mm	100 sheets
DRYMAT004	din	500 x 250mm	100 sheets
DRYMAT005	half din	250 x 250mm	100 sheets
DRYMAT006	din	500 x 280mm	100 sheets
DRYMAT007	double din	500 x 500mm	100 sheets
DRYMAT008	double din	600 x 600mm	100 sheets
DRYMAT009	double din	750 x 750mm	100 sheets
DRYMAT010	double din	900 x 900mm	100 sheets

GREEN DRY

Product code	Description	Size	Qty
B DRYMATGRN001	double din	750 x 500mm	50 sheets
DRYMATGRN002	double din	500 x 375mm	100 sheets
DRYMATGRN003	double din	610 x 305mm	100 sheets





Sterilisation Trays and Accessories - Protective Mats and Labels

GRID INSERTS AND MATS

Grid inserts and mats made of silicone. Heat resistant. Suitable for washer-disinfectors. Flexible.

Product code	Description	Colour	Heat-resistant	Size	Qty
A 08601	honeycomb structure, for blue perforated mesh trays	blue	up to +230°C	480 x 250mm	1 pc
B 08609	honeycomb structure, for blue perforated mesh trays	blue	up to +230°C	480 x 250mm	1 pc
C 08605	lattice structure	blue	-40°C up to +300°C	600 x 400 x 2mm	1 pc

SILICONE NUB MATS

Silicone nub mats. Suitable for washer-disinfectors. Flexible.

Product code	Description	Colour	Size	Qty
08426	compatible with ref 08422	blue	150 x 55mm	1 pc
08428	compatible with ref 08424	blue	200 x 135mm	1 pc
08429	compatible with ref 08425	blue	250 x 155mm	1 pc

SILICONE MESH COVER

Silicone mesh cover for system tray. Heat resistant up to 134°C. Flexible.

Product code	Description	Size	Qty
D 08881	for sieve trays 1/1	430 x 230mm/ 415 x 210mm	1 pc

STAINLESS STEEL SIEVE TRAY LABELS

Sieve tray labels are vital since hospitals are constantly in search of ever more sophisticated means of tracking sieve trays. Sieve tray labels are an ideal solution to this problem in terms of both colour coding and providing information. Heat resistant up to 134°C. For permanent use.

Product code	Description	Qty
E 01590	clip label, suitable for double frame	25 pcs
F 01445	label	50 pcs
G 01588	tie on label with hole (6 x 10mm)	5 pcs
H 01452	label with hole (4 x 9mm)	50 pcs



Sterilisation Trays and Accessories - Labels



LABEL HOLDER

Stainless Steel clip label holder for permanent use and heat resistant.

Product code	Description	Size	Qty
A 01497	label holder, compatible with ref 01550	50 x 20mm	25 pcs
B 01498	label holder, compatible with ref 01948	60 x 40mm	25 pcs

INSERTION LABEL

Compatible with clip label holder ref 01497 and dutiable label holders on containers. Heat resistant up to 134°C. For permanent use. Suitable for laser marking.

Product code	Colour	Size	Qty
C 01550	white	56 x 18mm	10 pcs
D 01551	red	56 x 18mm	10 pcs

ALUMINIUM LABEL

Suitable for slim label slots on containers and in ref 01497. Suitable for laser marking.

Product code	Colour	Size	Qty
E 28500	silver	18 x 50mm	25 pcs

