

## Comfort - Functionality - Durability

The Lumex Pivot-Arm Recliner features an easy-to-actuate Pivot-Arm to facilitate Recliner cleaning as well as patient entry and exit. The 22.5" wide seat and 400 lb weight capacity make the Lumex Pivot-Arm Recliner comfortable for a wider variety of patient types.

### Areas of Use:

- Dialysis
- Oncology
- IV Therapy
- PACU
- Emergency Room



Multiple Positions for Enhanced Patient Comfort

### Features:

- Both arms pivot outward, enabling both cleaning and easier patient transfers
- Three positions for patient comfort during treatment: Upright, First Recline and Second Recline
- Fourth position, caregiver-activated Trendelenburg Position (Shock Position), is foot activated via a red side pedal from either side of the Recliner
- The 22.5" wide seat ensures comfort for different patient types and sizes
- Weight capacity of 400 lb (181Kg), evenly distributed
- Four five-inch (5") Tente® Swivel Casters for ease of maneuverability in confined spaces. Three Total Locking Casters and one Directional Locking Caster
- Two self-storing side tables — with smooth surface and rounded corners for easy cleaning — provide patient prep and convenience
- Urethane armrest caps are comfortable for patient and easy to maintain; also available with fully upholstered arms
- Available in healthcare-grade vinyl or urethane upholstery in a broad palette of colors
- Weight-bearing, self-storing footplate for easier patient entry and exit and safer patient transport
- All Recliner parts are replaceable, greatly enhancing its lifespan and reducing its cost-in-use over time
- Meets:
  - California Technical Bulletin CA TB-117 Flammability Test
  - ANSI / BIFMA X5.4-2012
  - FNAE-80-214A, Section 7.8.1

### Dimensions

Overall Height..... 51.5" / 132 cm  
 Overall Width..... 34.5" / 88 cm  
 Overall Depth..... 38" / 97 cm  
 Depth Fully Reclined..... 68" / 173 cm  
 Seat Height ..... 22" / 56 cm  
 Seat Width ..... 22.5" / 57 cm  
 Seat Depth..... 19" / 48 cm  
 Back Height..... 32.5" / 83cm  
 Maximum Weight Capacity  
 Evenly Distributed..... 400 lb / 181 kg

### Options:

- Heat and Massage
- Fully Upholstered Armrests
- Headrest Skirt
- Headrest Cover
- Universal Pillow



Fully Weight-Bearing Footplate



Easy-to-operate yellow Pivot-Arm Release Lever on both sides of Recliner



Two Self-Storing Side Tables are standard