

MODEL : MET 6600 III

TRAUMA / TREATMENT / RECOVERY TROLLEY



| | |
|--------------------------|--|
| OVERALL DIMENSION | : L2140 x W720 x H600-900mm |
| MATERIAL | : MILDSTEEL FRAMEWORK |
| FINISHING | : EPOXY POWDER COATING |
| NETT WEIGHT | : 130 KG |
| LOAD CAPACITY | : < 280 KG |
| CASTOR DIAMETER | : Ø8 / 200mm CENTRAL LOCKING CUM Fifth Wheel |

STANDARD FEATURES :-

- Backrest elevated by gas lift (0 - 90°) and act as CPR release
- Hi-Lo elevated by hydraulic foot pedal
- Trendelenburg (18°) and Reverse Trendelenburg (-18°) adjustment by foot Pedal on both sides
- Foldable handles & Corner Bumper at head & foot end of trolley
- Bakelite mattress base for easy cleaning cum with restraining straps holders
- Full Length Sliding X-Ray translucent cassette holder
- 200mm swivel castors with central & directional locking with 5th wheel maneuver
- C arm compatible

INCLUSIVE OF :-

- 1 pair x Collapsible $\frac{3}{4}$ length chromed finish safety side rails with latch lock
- 1 x 75mm thick foam mattress with Cover
- 1 x Oxygen tank holder ('E' size)
- 2 x Mounted stainless steel IV pole with double hook
- 1 x Utility tray
- Drainage bag holder bars on both sides of bed