



neuINDUSTRIES
SAFER • CLEANER • BETTER

SEVERO

The best professional pill grinder for hospitals, nursing homes and home care organizations.



Made in The Netherlands

SEVERO Electric Pill Grinder

Key Features

- ✓ Fast and effective crushing of medication
- ✓ Avoiding exposure to medication and medication dust
- ✓ Small and portable design
- ✓ Easy to Clean
- ✓ Near silent operation
- ✓ Rechargeable battery
- ✓ Produced in the Netherlands

Products/Item codes

Severo Electric pill grinder (Item Code: EPC0031), including

- Stainless steel grinding head
- Power supply
- Optional: Stability Plate (STP0031)

Disposables

- Grinding cups (GRC2100): 2100 pcs/ctn
- PP Grinding Sheets(GRS2100): 2100 pcs/ctn
- Optional: Lids for Cup (GRL2000): 2000 pcs/ctn

Grinding Specifications

Grinding type	Autonomous grinding cycles via combined compression and rotation.
Pill grinding	Grinding of multiple standard pills and coated pills per cycle
Grinding Cup	30 ml (Ø45 x H28 mm) non-grooved inner surface prevents powder entrapment.
Grinding cycle time	8 - 14 seconds
Sound level	52 dB (soft whisper level)
Battery last	1000 cycles
Charging time	2 hours
Switches	On/Off
Indicators	Battery load and battery quality indicator

Technical Specifications

Nominal Voltage	12 Volts
Approximate weight	1.5 Kilograms
Dimensions (LxWxH)	117 x 86 x 232 mm
Terminals	5,5 x 2.1 x 12 mm DC Connector for recharging
AC/DC adaptor	100-240 VAC input, 12 VDC output
Battery	Rechargeable sealed lithium-ion, 12 Volts, 3 Ah
Battery self-discharge	1% of capacity per month (25 c°)

Quality Assurance

Before shipment each device is thoroughly tested/inspected on:

- ✓ Scratches, dents, dirt or other exterior deformations
- ✓ Odd sounds during operation
- ✓ Functioning and grinding quality
- ✓ Battery capacity
- ✓ Availability of accessories

For further information



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.neuindustries.com | sales@neuindustries.com



IMS Medical B.V
www.severo.nl