


Surgical Tapes

Secure. Gentle. Reliable



Silicone Tape

Silicone Tape features a gentle, hypoallergenic silicone adhesive that securely holds dressings, fix tubes, catheters, and medical devices without irritating the skin. Its breathable and waterproof design protects wounds from contamination while allowing the skin to breathe and heal. Ideal for sensitive and delicate skin, it provides comfortable fixation with easy, painless removal, leaving no residue. This tape is perfect for surgical, nursing, and home care applications.



Specifications

Material: Silicone

Color: Transparent






Shelf life: 3 years

Store in a cool, dry place. Keep away from direct sunlight.

Features

- Waterproof
- Breathable
- Painless removal with no residue
- Gentle silicone ideal for sensitive skin

User guide


- ①  Clean and dry the skin
- ②  Cut tape to desired length
- ③  Place tape over skin and dressing/ tubing
- ④  Secure IV Line using silicone tape
- ⑤  Secure Nasal Gastric Tube using silicone tape

Ordering Information

Item Code	Size (W x L)	Box Size	Carton Size	Packing
SIT001A	1.9cm x 1.5m	17 x 11.5 x 6.5cm	36 x 25 x 34.5cm	12 roll/box, 20 box/ctn
SIT001B	1.9cm x 5m	17 x 11.5 x 6.5cm	36 x 25 x 34.5cm	12 roll/box, 20 box/ctn
SIT002A	2.5cm x 1.5m	17 x 11.5 x 6.5cm	36 x 25 x 34.5cm	12 roll/box, 20 box/ctn
SIT002B	2.5cm x 5m	17 x 11.5 x 6.5cm	36 x 25 x 34.5cm	12 roll/box, 20 box/ctn
SIT003A	5cm x 1.5m	17 x 11.5 x 6cm	36 x 25 x 32cm	6 roll/box, 20 box/ctn
SIT003B	5cm x 5m	17 x 11.5 x 6cm	36 x 25 x 32cm	6 roll/box, 20 box/ctn
SIT004A	7.5cm x 1.5m	17 x 11.5 x 8.5cm	36 x 25 x 19cm	6 roll/box, 8 box/ctn
SIT004B	7.5cm x 5m	17 x 11.5 x 8.5cm	36 x 25 x 19cm	6 roll/box, 8 box/ctn

PU Fixing Tape

PU Fixing Tape is a transparent, waterproof, and breathable medical adhesive film designed for secure fixation of dressings, catheters, and medical devices. Made from polyurethane (PU) film with a skin-friendly adhesive, it conforms smoothly to the skin, allowing natural movement while protecting the site from water, bacteria, and friction.



Specifications

Material: Polyurethane fabric

Color: Transparent


Shelf life: 3 years


Store in a cool, dry place. Keep away from direct sunlight.


Features


- Waterproof & shower-safe
- Bacteria-resistant
- Prevents friction & leakage
- Breathable
- Skin-safe & painless removal


User guide

- 

1 Cut the tape to desired length
- 

2 Peel off one side of the protective film
- 

3 Apply tape and peel off the remaining film
- 

4 Flip and peel off the remaining film
- 

5 Apply tape and press gently

Ordering Information

Item Code	Size (W x L)	Box Size	Carton Size	Packing
PUT001	5.0cm x 10m	8 x 8 x 5.5cm	42 x 42 x 35cm	1roll/box, 150 box/ctn
PUT002	7.5cm x 10m	8 x 8 x 7.8cm	42 x 42 x 41cm	1roll/box, 125 box/ctn
PUT003	10cm x 10m	8 x 8 x 10.5cm	42 x 34 x 44cm	1roll/box, 80 box/ctn
PUT004	15cm x 10m	8 x 8 x 15.5cm	42 x 42 x 33cm	1roll/box, 50 box/ctn
PUT005	20cm x 10m	8 x 8 x 20.5cm	42 x 42 x 43cm	1roll/box, 50 box/ctn
PUT006	30cm x 10m	8 x 8 x 30.5cm	42 x 42 x 32.5cm	1roll/box, 25 box/ctn

Non-woven Fixing Tape

Non-Woven Fixing Tape is a soft, breathable adhesive tape designed to securely fix medical devices such as dressings, syringes, catheters, and gauze. Made from non-woven fabric, it is lightweight, flexible, and gentle on the skin, making it ideal for everyday clinical and nursing applications. The tape allows moisture and air to pass through while providing reliable fixation without irritating the skin.



Specifications

Material: Non-woven fabric

Color: White


Shelf life: 3 years

Store in a cool, dry place. Keep away from direct sunlight.

Features

- Soft and skin-friendly material
- Easy to cut and apply
- Breathable

User guide



1 Place primary dressing on the required area

2 Cut tape to desired length

3 Peel off one side of the protective film

4 Apply tape and peel off the remaining film

5 Press gently to secure

Ordering Information

Item Code	Size (W x L)	Box Size	Carton Size	Packing
FIXT001	5.0cm x 10m	8 x 8 x 5.5cm	42 x 42 x 35cm	1 roll/box, 150 box/ctn
FIXT002	7.5cm x 10m	8 x 8 x 7.8cm	42 x 42 x 41cm	1 roll/box, 125 box/ctn
FIXT003	10cm x 10m	8 x 8 x 10.5cm	42 x 34 x 44cm	1 roll/box, 80 box/ctn
FIXT004	15cm x 10m	8 x 8 x 15.5cm	42 x 42 x 33cm	1 roll/box, 50 box/ctn
FIXT005	20cm x 10m	8 x 8 x 20.5cm	42 x 42 x 43cm	1 roll/box, 50 box/ctn
FIXT006	30cm x 10m	8 x 8 x 30.5cm	42 x 42 x 32.5cm	1 roll/box, 25 box/ctn

Rev1 (0825)

