

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-AUB1000

Description	Adaptor Urinal Bulbous Disposable Pulp	
Ref. No.	AUB1000	
UDI-DI	95510204301752	
Product Image		
Rated Capacity	Not applicable	
Nominal Weight	9g +/- 2g	
Dimension	145 x 115 x 65mm	
Graduated Marking	Not applicable	
Material	Biodegradable Kraft	
Colour	Golden Brown	
Base	None	
Finishing	Matt	
Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"	
Fluid Temperature	Initial 35°C	
Fluid Retention	> 4 hours	
Packing	200 pieces per carton	
Carton Weight	2.5 kg	
Carton Size	L510 x W260 x H250 mm	
Disposal	Macerator or appropriate waste bin	
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.	
Shelf Life	7 years	

Associated Accessories

Moulded Pulp	MUB1000	Male Urinal Bulbous
--------------	---------	---------------------

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com