

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-AUB1000

Description	Adaptor Urinal Bulbous Disposable Pulp
Ref. No.	AUB1000
UDI-DI	95510204301752
Product Image	
Rated Capacity	Not applicable
Nominal Weight	9g +/- 2g
Dimension	145 x 115 x 65mm
Graduated Marking	Not applicable
Material	Biodegradable Kraft
Colour	Golden Brown
Base	None
Finishing	Matt
Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Fluid Retention	> 4 hours
Packing	200 pieces per carton
Carton Weight	2.5 kg
Carton Size	L510 x W260 x H250 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Associated Accessories

Moulded Pulp	MUB1000	Male Urinal Bulbous
--------------	---------	---------------------

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

