


TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-BPG2000

Description	Bedpan Liner Graduated 2000ml Disposable Pulp	
Ref. No.	BPG2000	
UDI-DI	95510204300076	
Product Image		
Rated Capacity	2000 ml	
Nominal Weight	64g +/- 5g	
Dimension	360 x 300 x 100 mm	
Graduated Marking	Minimum 100 ml ; Maximum 2000 ml ; Interval 200ml	
Material	Biodegradable Kraft	
Colour	Golden Brown	
Base	Flat, Stable	
Finishing	Matt	
Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"	
Fluid Temperature	Initial 35°C	
Fluid Retention	> 4 hours	
Packing	100 pieces per carton	
Carton Weight	7.5 kg	
Carton Size	L755 x W370 x H310 mm	
Disposal	Macerator or appropriate waste bin	
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.	
Shelf Life	7 years	

Associated Accessories

Moulded Pulp	BPS2000	Bedpan Support < 100 kg
Plastic Support	SBP2000	Antimicrobial Bedpan Support < 160 kg

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com