

## TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-CUP0090

Description Gallipot Food Grade 90ml Disposable Pulp  
Ref. No. CUP0090  
UDI-DI 95510204302506  
Product Image



Rated Capacity 90 ml  
Nominal Weight 5.3g +/- 1g  
Dimension Ø75 x 50 mm  
Graduated Marking None

Material Biodegradable Sugarcane Bagasse  
Colour White  
Base Flat, Stable  
Finishing Smooth internal ; Mesh external

Manufacturing Compliance FDA 21 CFR 176.170  
EU Regulation No. 1935/2004  
PAS29:1999 "Disposable Pulp Products For Use in Healthcare"

Fluid Temperature Initial 35°C  
Fluid Retention > 4 hours

Packing 1000 pieces per carton  
Carton Weight 5.8 kg  
Carton Size L480 x W400 x H390 mm

Disposal Macerator or appropriate waste bin  
Storage Condition Store in cool, dry environment away from open flames and other ignition source.  
No direct sunlight.  
Shelf Life 7 years

**Disclaimer:**

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



**NEU Industries Pte Ltd**

132 Joo Seng Rd, #04-01, Singapore 368358  
Tel: +65 6817 8969 Fax: +65 6636 9813  
www.neuindustries.com | sales@neuindustries.com