

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-CUP0300

Description	Specimen Cup 300ml Disposable Pulp
Ref. No.	CUP0300
UDI-DI	95510204302438
Product Image	
Rated Capacity	300 ml
Nominal Weight	10g +/- 1g
Dimension	Ø93 x 92 mm
Graduated Marking	None
Material	Biodegradable Sugarcane Bagasse
Colour	White
Base	Flat, Stable
Finishing	Smooth internal ; Mesh external
Manufacturing Compliance	FDA 21 CFR 176.170 EU Regulation No. 1935/2004 PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Fluid Retention	> 4 hours
Packing	1000 pieces per carton
Carton Weight	10.5 kg
Carton Size	L480 x W450 x H395 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

