

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-GPB1000

Description General Purpose Bowl 1000ml Disposable Pulp
Ref. No. GPB1000
UDI-DI 95510204300694
Product Image



Rated Capacity 1000 ml
Nominal Weight 29g +/- 5g
Dimension 220 x 170 x 85 mm
Graduated Marking Minimum 200ml ; Maximum 1000ml ; Interval 200ml

Material Biodegradable Kraft
Colour Golden Brown
Base Flat, Stable
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 100 pieces per carton
Carton Weight 3.5 kg
Carton Size L510 x W260 x H360 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

