


TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-GPB3000

Description	General Purpose Bowl 3000ml Disposable Pulp
Ref. No.	GPB3000
UDI-DI	95510204300830
Product Image	
Rated Capacity	3000 ml
Nominal Weight	64g +/- 5g
Dimension	300 x 300 x 100 mm
Graduated Marking	Minimum 500ml ; Maximum 3000ml ; Interval 500ml
Material	Biodegradable Kraft
Colour	Golden Brown
Base	Flat, Stable
Finishing	Matt
Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Fluid Retention	> 4 hours
Packing	100 pieces per carton
Carton Weight	7.0 kg
Carton Size	L755 x W370 x H310 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev2(0825)

