

## TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-GPP2000

Description	General Purpose Pot 2000ml Disposable Pulp
Ref. No.	GPP2000
UDI-DI	95510204301134
Product Image	
Rated Capacity	2000 ml
Nominal Weight	42g +/- 5g
Dimension	260 x 180 x 105 mm
Graduated Marking	Minimum 50ml ; Maximum 2000ml ; Interval 100ml
Material	Biodegradable Kraft
Colour	Golden Brown
Base	Flat, Stable
Finishing	Matt
Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Fluid Retention	> 4 hours
Packing	100 pieces per carton
Carton Weight	5.5 kg
Carton Size	L755 x W370 x H310 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

### Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



**NEU Industries Pte Ltd**

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com