

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-GSB1500

Description General Semi Bowl 1500ml Disposable Pulp
Ref. No. GSB1500
UDI-DI 95510204300762
Product Image



Rated Capacity 1500 ml
Nominal Weight 45g +/- 5g
Dimension 330 x 220 x 110 mm
Graduated Marking Minimum 50ml ; Maximum 1400ml ; Interval 100ml

Material Biodegradable Kraft
Colour Golden Brown
Base Flat, Stable
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 100 pieces per carton
Carton Weight 5.5 kg
Carton Size L450 x W340 x H450 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

