

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-GWB3000

Description General Wash Bowl 4000ml Disposable Pulp
Ref. No. GWB3000
UDI-DI 95510204302810
Product Image



Rated Capacity 3000 ml
Nominal Weight 59g +/- 5g
Dimension 300 x 300 x 100 mm
Graduated Marking Minimum 500ml ; Maximum 3000ml ; Interval 500ml

Material Biodegradable Kraft
Colour Golden Brown
Base Flat, Stable
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature 45°C +/- 3°C
Fluid Retention > 4 hours (depending on type of detergent used)



Packing 100 pieces per carton
Carton Weight 7.5 kg
Carton Size L755 x W370 x H310 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

