

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-GWB4000

Description General Wash Bowl 4000ml Disposable Pulp
Ref. No. GWB4000
UDI-DI 95510204300908
Product Image



Rated Capacity 4000 ml
Nominal Weight 68g +/- 5g
Dimension 350 x 255 x 120 mm
Graduated Marking Minimum 500ml ; Maximum 4000ml ; Interval 500ml

Material Biodegradable Kraft
Colour Golden Brown
Base Flat, Stable
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours
Detergent-Proof Max. 2 hours, water with soap at 45°C



Packing 100 pieces per carton
Carton Weight 7.5 kg
Carton Size L520 x W360 x H500 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358
Tel: +65 6817 8969 Fax: +65 6636 9813
www.neuindustries.com | sales@neuindustries.com