

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-KND0700

Description Kidney Dish 700ml Disposable Pulp
Ref. No. KND0700
UDI-DI 95510204301370
Product Image



Rated Capacity 700 ml
Nominal Weight 24g +/- 5g
Dimension 270 x 140 x 50 mm
Graduated Marking None

Material Biodegradable Kraft
Colour Golden Brown
Base Flat, Stable
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 300 pieces per carton
Carton Weight 7.5 kg
Carton Size L465 x W275 x H320 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com