

## TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-KND0950

Description Kidney Dish 950ml Disposable Pulp  
Ref. No. KND0950  
UDI-DI 95510204301448  
Product Image



Rated Capacity 950 ml  
Nominal Weight 29g +/- 5g  
Dimension 270 x 140 x 55 mm  
Graduated Marking Minimum 100 ml ; Maximum 900 ml ; Intervals of 200 ml

Material Biodegradable Kraft  
Colour Golden Brown  
Base Flat, Stable  
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"  
Fluid Temperature Initial 35°C  
Fluid Retention > 4 hours

Packing 300 pieces per carton  
Carton Weight 8.6 kg  
Carton Size L580 x W360 x H420 mm

Disposal Macerator or appropriate waste bin  
Storage Condition Store in cool, dry environment away from open flames and other ignition source.  
No direct sunlight.  
Shelf Life 7 years

### Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



**NEU Industries Pte Ltd**

132 Joo Seng Rd, #04-01, Singapore 368358  
Tel: +65 6817 8969 Fax: +65 6636 9813  
www.neuindustries.com | sales@neuindustries.com