

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-KND0950

Description Kidney Dish 950ml Disposable Pulp
Ref. No. KND0950
UDI-DI 95510204301448
Product Image



Rated Capacity 950 ml
Nominal Weight 29g +/- 5g
Dimension 270 x 140 x 55 mm
Graduated Marking Minimum 100 ml ; Maximum 900 ml ; Intervals of 200 ml

Material Biodegradable Kraft
Colour Golden Brown
Base Flat, Stable
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 300 pieces per carton
Carton Weight 8.6 kg
Carton Size L580 x W360 x H420 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

