

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-KNS0950

Description Kidney Dish Food Grade 950ml Disposable Pulp
Ref. No. KNS0950
UDI-DI 95510204301516
Product Image



Rated Capacity 950 ml
Nominal Weight 23g +/- 2g
Dimension 245 x 120 x 45 mm
Graduated Marking None

Material Biodegradable Sugarcane Bagasse
Colour White
Base Flat, Stable
Finishing Smooth internal ; Mesh external

Manufacturing Compliance FDA 21 CFR 176.170
EU Regulation No. 1935/2004
PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 300 pieces per carton
Carton Weight 7.4 kg
Carton Size L470 x W275 x H430 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358
Tel: +65 6817 8969 Fax: +65 6636 9813
www.neuindustries.com | sales@neuindustries.com