


TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-LBP2000

Description	Lid for Bedpan Liner and Midi Slipper Pan Liner
Ref. No.	LBP2000
UDI-DI	95510204302988
Product Image	
Nominal Weight	22g +/- 1g
Dimension	330 x 287 x 0.4 mm
Material	Paper
Colour	White
Finishing	Matt
Manufacturing Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Packing	200 pieces per carton
Carton Weight	4.6 kg
Carton Size	L340 x W300 x H90 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com