

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-LGP2000

Description
Ref. No.
UDI-DI
Product Image

Lid for General Purpose Pot
LGP2000
95510204303046



Nominal Weight
Dimension

8.3g +/- 1g
276 x 146 x 0.4 mm

Material
Colour
Finishing

Paper
White
Matte

Manufacturing Compliance
Fluid Temperature

PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Initial 35°C

Packing
Carton Weight
Carton Size

200 pieces per carton
1.7 kg
L290 x W160 x H90 mm

Disposal
Storage Condition

Macerator or appropriate waste bin
Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.

Shelf Life

7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)

