

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-LGP2000

Description	Lid for General Purpose Pot
Ref. No.	LGP2000
UDI-DI	95510204303046
Product Image	
Nominal Weight	8.3g +/- 1g
Dimension	276 x 146 x 0.4 mm
Material	Paper
Colour	White
Finishing	Matte
Manufacturing Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Packing	200 pieces per carton
Carton Weight	1.7 kg
Carton Size	L290 x W160 x H90 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

