

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-LMU1000

Description Lid for Male Urinal Bulbous Disposable Pulp
Ref. No. LMU1000
UDI-DI 95510204303114
Product Image



Rated Capacity Not applicable
Nominal Weight 4g +/- 1g
Dimension 85 x 70 x 15 mm
Graduated Marking Not applicable

Material Biodegradable Kraft
Colour Golden Brown
Base Flat, Stable
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 1000 pieces per carton
Carton Weight 4.0 kg
Carton Size L510 x W260 x H360 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
Shelf Life No direct sunlight.
7 years

Associated Accessories

Moulded Pulp MUB1000 Male Urinal Bulbous

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358
Tel: +65 6817 8969 Fax: +65 6636 9813
www.neuindustries.com | sales@neuindustries.com