

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-MSP1300

Description Midi Slipper Pan Liner 1300ml Disposable Pulp
Ref. No. MSP1300
UDI-DI 95510204300212
Product Image



Rated Capacity 1300 ml
Nominal Weight 58g +/- 5g
Dimension 345 x 305 x 80 mm
Graduated Marking Minimum 200 ml ; Maximum 1200 ml ; Interval 200 ml

Material Biodegradable Kraft
Colour Golden Brown
Base Flat, Stable
Finishing Matt

Compliance PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 100 pieces per carton
Carton Weight 7.0 kg
Carton Size L580 x W360 x H420 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
Shelf Life No direct sunlight.
7 years

Associated Accessories

Plastic Support SSP1300 Antimicrobial Midi Slipper Pan Support < 160 kg

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

