

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-MSP1300

Description	Midi Slipper Pan Liner 1300ml Disposable Pulp
Ref. No.	MSP1300
UDI-DI	95510204300212
Product Image	
Rated Capacity	1300 ml
Nominal Weight	58g +/- 5g
Dimension	345 x 305 x 80 mm
Graduated Marking	Minimum 200 ml ; Maximum 1200 ml ; Interval 200 ml
Material	Biodegradable Kraft
Colour	Golden Brown
Base	Flat, Stable
Finishing	Matt
Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Fluid Retention	> 4 hours
Packing	100 pieces per carton
Carton Weight	7.0 kg
Carton Size	L580 x W360 x H420 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Associated Accessories

Plastic Support	SSP1300	Antimicrobial Midi Slipper Pan Support < 160 kg
-----------------	---------	---

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com