

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-MSP2000

Description	Maxi Slipper Pan Liner 2000ml Disposable Pulp
Ref. No.	MSP2000
UDI-DI	95510204300458
Product Image	
Rated Capacity	2000 ml
Nominal Weight	69g +/- 5g
Dimension	375 x 275 x 95 mm
Graduated Marking	Minimum 200ml ; Maximum 2000ml ; Interval 400ml
Material	Biodegradable Kraft
Colour	Golden Brown
Base	Flat, Stable
Finishing	Matt
Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Fluid Retention	> 4 hours
Packing	100 pieces per carton
Carton Weight	7.5 kg
Carton Size	L590 x W395 x H560 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Associated Accessories

Plastic Support	SSP2000	Maxi Slipper Pan Support < 160 kg
-----------------	---------	-----------------------------------

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

