


TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-MUB1000

Description	Male Urinal Bulbous 1000ml Disposable Pulp	
Ref. No.	MUB1000	
UDI-DI	95510204301684	
Product Image		
Rated Capacity	1000 ml	
Nominal Weight	53g +/- 5g	
Dimension	245 x 110 x 125 mm	
Graduated Marking	None	
Material	Biodegradable Kraft	
Colour	Golden Brown	
Base	Flat, Stable	
Finishing	Matt	
Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"	
Fluid Temperature	Initial 35°C	
Fluid Retention	> 4 hours	
Packing	100 pieces per carton	
Carton Weight	7.0 kg	
Carton Size	L710 x W410 x H790 mm	
Disposal	Macerator or appropriate waste bin	
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.	

Associated Accessories

Moulded Pulp	LMU1000	Lid for Male Urinal Bulbous
	AUB1000	Adaptor Urinal Bulbous

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com