

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-T2C2000

Description	Anaesthetic Tray 2-Compartment 2000ml Disposable Pulp
Ref. No.	T2C2000
UDI-DI	95510204302056
Product Image	
Rated Capacity	2000 ml
Nominal Weight	28g +/- 2g
Dimension	225 x 180 x 50 mm
Graduated Marking	None
Material	Biodegradable Sugarcane Bagasse
Colour	White
Base	Flat, Stable
Finishing	Smooth internal ; Mesh external
Manufacturing Compliance	FDA 21 CFR 176.170 EU Regulation No. 1935/2004 PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Fluid Retention	> 4 hours
Packing	300 pieces per carton
Carton Weight	8.9 kg
Carton Size	L510 x W260 x H210 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com