

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-T4C2000

Description Anaesthetic Tray 4-Compartment 2000ml Disposable Pulp
Ref. No. T4C2000
UDI-DI 95510204301998
Product Image



Rated Capacity 2000 ml
Nominal Weight 35g +/- 2g
Dimension 270 x 180 x 45 mm
Graduated Marking None

Material Biodegradable Sugarcane Bagasse
Colour White
Base Flat, Stable
Finishing Smooth internal ; Mesh external

Manufacturing Compliance FDA 21 CFR 176.170
EU Regulation No. 1935/2004
PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 100 pieces per carton
Carton Weight 4.0 kg
Carton Size L370 x W285 x H185 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358
Tel: +65 6817 8969 Fax: +65 6636 9813
www.neuindustries.com | sales@neuindustries.com