

TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-TIS1000

Description Instrument Tray 1000ml Disposable Pulp
Ref. No. TIS1000
UDI-DI 95510204302124
Product Image



Rated Capacity 1000 ml
Nominal Weight 21g +/- 2g
Dimension 195 x 135 x 48 mm
Graduated Marking None

Material Biodegradable Sugarcane Bagasse
Colour White
Base Flat, Stable
Finishing Smooth internal ; Mesh external

Manufacturing Compliance FDA 21 CFR 176.170
EU Regulation No. 1935/2004
PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature Initial 35°C
Fluid Retention > 4 hours

Packing 500 pieces per carton
Carton Weight 11.0 kg
Carton Size L420 x W210 x H330 mm

Disposal Macerator or appropriate waste bin
Storage Condition Store in cool, dry environment away from open flames and other ignition source.
No direct sunlight.
Shelf Life 7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

