


TECHNICAL DATA SPECIFICATIONS

Doc. No. NEU-TDS-UMT1000

Description	Urinal Measuring Tape
Ref. No.	UMT1000
UDI-DI	Pending
Product Image	
Nominal Weight	2.3g +/- 0.1g
Dimension	300 x 15 x 0.4 mm
Graduated Marking	Minimum 100 ml ; Maximum 1000 ml ; Intervals of 100 ml
Material	Paper
Colour	Light Brown
Base	Flat, Stable
Finishing	Matte
Manufacturing Compliance	PAS29:1999 "Disposable Pulp Products For Use in Healthcare"
Fluid Temperature	Initial 35°C
Packing	1000 pieces per pack
Pack Weight	2.4 kg
Pack Size	L300 x W150 x H70 mm
Disposal	Macerator or appropriate waste bin
Storage Condition	Store in cool, dry environment away from open flames and other ignition source. No direct sunlight.
Shelf Life	7 years

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Technical Data Specifications or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0725)



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com