

Pink Card Bowie & Dick Test



The Pink Card is a Bowie & Dick test designed to detect air leaks, inadequate steam penetration and vacuum pump failures in pre-vacuum steam sterilizers. The test is conducted each day in an empty sterilizer before processing any loads.

- ★ Water-based and non-toxic chemical indicator is used.
- ★ It is easy to read and interpret.
- ★ The steel apparatus (BD-8 HLD) of the Pink Card allows the Pink Card to hold the card in place during the sterilization cycle.
- ★ The metal ring apparatus (BD-8 RNG) of the Pink Card allows Pink Card to be attached to a rack in the sterilizer and holds the card in place during the sterilization cycle.
- ★ Compatible with ISO 11140-1 and EN 867-5.
- ★ Different designs are made for the company/person on request.


BD-8 HLD



BD-8 PNG



BD-8 FRONT SIDE



BD-8 BACK SIDE



PRODUCT SIZE: 50 X 90 mm

REF NO	DESCRIPTION	PACKAGING
BD-8	Pink Card Bowie & Dick Test (Pack of 100)	60 Packs / CARTON
BD-8 HLD	Holder for BD-8 Pink Card (Stainless Steel)	45 Pcs. / CARTON
BD-8 RNG	Metal Ring for BD-8 Pink Card	100 Pcs. / CARTON