

ELECTRIC TRANSFER SHOWER COMMUNE CHAIR



Introducing the NEW revolutionary ELECTRIC Transfer Shower Commode Chair that will be a great relieve to the Caregivers!

Multiple Functions in One chair

Transfer

Transfer residents between bed, sofa, dining table and toilet bowl effortlessly!

Shower

Water resistance IP65 rating for showering. Removable controller to prevent water ingress.

Commode

Fits over all toilet bowls with electric seat height up and down adjustment at your finger tips.

Transfer from Bed to Toilet Bowl



| Description | Specification |
|-------------------|------------------------|
| Material | Robust steel structure |
| Lift & Descend | Electric |
| Seat Depth | 395mm |
| Seat Width | 485mm |
| Seat Height | 410 to 605 mm |
| Seat Opening | 180° outwards |
| Total Length | 830mm |
| Total Width | 525mm |
| Total Height | 835 to 1035 mm |
| Gap Clearance | 120mm |
| Castor wheels | Front 125mm, Rear 75mm |
| Electric supply | 24VDC @ 5A |
| Power | 120W |
| Battery | 2500mAh |
| Net weight | 28.5kg |
| Safe Working Load | 100kg |
| Water Resistance | IP65 |



+65 9820 9495



sales@neuindustries.com



www.neuindustries.com

NEU Industries Pte Ltd

132 Joo Seng Road #04-01. Singapore 368358

T: +65 6817 8969 | F: +65 6636 9813

