

MANUAL CRANK TRANSFER SHOWER COMMODE CHAIR

ESSENTIAL
AID FOR
CAREGIVERS



MANUAL CRANK TRANSFER SHOWER COMMODE CHAIR

Model: SCM2000

Shower Commode Manual

Introducing the NEW revolutionary MANUAL Crank Transfer Shower Commode Chair that will be a great relieve to the Caregivers!

Multiple Functions in One Chair

Transfer

Transfer residents between bed, sofa, dining table and toilet bowl effortlessly!

Shower

Water resistance IP65 rating for showering. Removable controller to prevent water ingress.

Commode

Fits over all toilet bowls with electric seat height up and down adjustment at your finger tips.

Transfer from Bed to Toilet Bowl



Description	Specification
Material	Robust steel structure
Lift & Descend	Manual Crank
Seat Depth	420mm
Seat Width	465mm
Seat Height	450 to 650 mm
Seat Opening	180° outwards
Total Length	700mm
Total Width	595mm
Total Height	805 to 1005 mm
Gap Clearance	148mm
Castor wheels	100mm
Net weight	25.5kg
Safe Working Load	100kg
Water Resistance	IP65



+65 9820 9495



sales@neuindustries.com



www.neuindustries.com

NEU Industries Pte Ltd

132 Joo Seng Road #04-01. Singapore 368358

T: +65 6817 8969 | F: +65 6636 9813

