

MODEL : SCS2000 Series

STANDARD SHOWER cum COMMODE MOBILE CHAIR



Ordering Information

SCS2000-NOB-CH

SCS2000-FRB-CH

SCS2000-REB-CH

NOB = **NO** SAFETY SEAT **BELT**

FRB = **FRONT BUCKLE** SAFETY SEAT **BELT**

REB = **REAR BUCKLE** SAFETY SEAT **BELT**

CH = **CLEARNCE HEIGHT** OVER TOILET BOWL

Standard Item

Plastic Commode Pan with Cover



Optional Item

Stainless Steel Holder (Removable) under seat for Disposable Bedpan Liner



OVERALL DIMENSION	: W510 x D630/850 (footrest) x H1000mm
SEAT DIMENSION	: W450 x D460mm
CLEARANCE HEIGHT (CH)	: To specify height over toilet bowl upon ordering
MATERIAL	: Rust and Corrosion Resistant Stainless Steel Frame
WEIGHT CAPACITY	: Maximum 150kg

STANDARD FEATURES

- Stainless steel frame
- Smooth push handle for easy transportation and cleaning
- Removable armrest with anti-slip surface
- Retractable footrest with anti-tipper
- ABS Plastic on backrest, seat & armrest
- Easy to wheel over standard toilet bowl as elevated toilet seat
- 4 x 125 mm swivel castor with individual brake

OPTIONS

I.V. Drip Pole	: Removable fixed height with double hook (Model: SCSIVP2)
Oxygen Tank Holder	: Removable mounted at rear (Model: SCSOTH)