

MODEL : SCS2000 Series

STANDARD SHOWER cum COMMUNE MOBILE CHAIR



Ordering Information

SCS2000**BPL**-Options

SCS2000**PCP**-Options

BPL

Stainless Steel Holder (Removable)
for Disposable Bedpan Liner (**BPL**)



PCP

Plastic Commode Pan (**PCP**)



Options

Extended Seat Height (**-EH**)

Safety Seat Belt Front Buckle (**-FB**)

Safety Seat Belt Rear Buckle (**-RB**)

OVERALL DIMENSION

: W510 x D630/850 (footrest) x H1000mm

SEAT DIMENSION / HEIGHT

: W450 x D460mm / H545mm (Std) or H565mm (Ext)

CLEARANCE OVER WC

: H460mm (Standard) / H480mm (Extended)

Other customized options available

MATERIAL

: Rust and Corrosion Resistant Stainless Steel Frame :

WEIGHT CAPACITY

Maximum 150kg

STANDARD FEATURES

- Smooth push handle for easy transportation and cleaning
- Removable armrest with Anti-Slip surface
- Retractable footrest with Anti-Tipper
- Stainless Steel holder for commode pan, bedpan or pulp
- ABS Plastic on backrest, seat & armrest
- Easy to wheel over standard toilet bowl as elevated toilet seat
- 4 x 125 mm Swivel Castor with Brake
- Stainless Steel Frame

OPTIONS

Safety Seat Belt : Front buckle (-FB) or rear buckle (-RB)

I.V. Drip Pole : Removable fixed height with double hook (Model: SCSIVP2)

Height Over WC : Extended height (-EH) over WC customizable