

ELECTRIC DIALYSIS CHAIR


SKE-133



PRODUCT DATA SHEET



SKE-133



Technical Paramters

▪ External Size (LxW):	2000 x 580 mm (±10 mm)
▪ Back:	770 x 580 mm (±10 mm)
▪ Seat:	470 x 580 mm (±10 mm)
▪ Leg:	400 x 540 mm (±10 mm)
▪ Foot Board:	430 x 210 mm (±10 mm)
▪ Cushion Thickness:	80 mm (±10 mm)
▪ Safe Working Load:	200 kg

Functions

▪ Backrest Adjustment:	-8° - 70° (±5°)
▪ Leg Rest Adjustment:	-60° - 0° (±5°)
▪ Height Adjustment:	550 - 700 mm (±10 mm)
▪ With CPR	
▪ Armrest Rotation:	180° Horizontal Swing

Technical Configuration

☑ Hand Controller:	1 pc
☑ 3-inch Dual-side Castors:	4 pcs, all with brakes
☑ Armrest:	2 pcs
☑ Side Table:	1 pc (shown in picture)



FEATURES



Steel
Baosteel from Fortune Global 500.



Plastic
11-stage epoxy coating process with a 0.12 mm paint thickness, ASTM-tested antibacterial properties, and impact resistance of up to 50 kg.




Welding
Panasonic Robotic ensure 360° full smooth welding.



Electric Motor
Linear motor compliant with UL, RoHS, and EN standards, operating below 50 dB with a rated lifespan of 20,000 cycles.




Chair Upholstery
Available in genuine leather or supebric upholstery, featuring quality craftsmanship and a soft, comfortable finish for enhanced patient comfort.



Foot Rest
The footrest automatically aligns with the raised legrest, creating a flat support surface and extending the patient area by 250 mm.



Ergonomic Armrest
The pivoting armrest allows flexible positioning, while its concave design provides stable support, prevents arm slippage, and enhances patient comfort.



ABS Tray Table
Attachable tray table with anti-skid convex edges. Has built-in cup holder and phone holder for patient convenience.

Ultrafiber Leather Colour Table



* Note: the above is a physical scan and there may be certain color differences. Please refer to the actual product for accurate colors.



NEU Industries Pte Ltd

132 Joo Seng Rd, #04-01, Singapore 368358

Tel: +65 6817 8969 | Fax: +65 6636 9813

www.neuindustries.com | sales@neuindustries.com

