

Simply AutoFit™



World's first AutoFit mask

F&P Solo™ combines years of relentless innovation and cutting-edge components to usher in a new era of automatic masks.

The world's first AutoFit mask, F&P Solo, simplifies mask setup like no other, using AutoLock™ Technology.



Fit

Users simply stretch the mask over their head and onto their face.



Touch

With one touch, users adjust for their preferred fit by sealing the mask under their nose with a gentle push.



AutoLock technology holds the seal in place and requires only a single pull to release the headgear.

Sleep

AutoLock technology comfortably holds the mask in place, dynamically responding to sleep movement and CPAP-therapy pressures.



One revolutionary headgear design

Simple,
every night



Users can easily fine-tune the mask during nightly use for an optimal fit.

X-ray view

Reduces
overtightening




Without conventional straps to undo and pull, the AutoFit headgear reduces overtightening.

Quality
components




AutoLock technology,
manufactured with carbon fibre
material for seamless movement.

Two comfortable cushion options



Nasal Cushion

Sits under the nose allowing for freedom of movement while keeping the cushion comfortably in place.



Pillows Cushion

Soft, thin pillows cushion sits comfortably inside your nose providing a stable seal.



Simplify the mask experience across the patient's journey



Mask Selection

One headgear system with two cushion options to simplify mask choice for you and your patient.




Setup with the F&P myMask™ App

For an enhanced experience, F&P myMask supports your patients with mask setup and fitting.




User-friendly packaging

Step-by-step pictorial instructions provides confidence to set up the mask independently.



F&P myMask support

F&P myMask helps with fitting, cleaning and troubleshooting. Support when it's needed most.



Simple for you and your patients

F&P Solo has done the hard work so you and your patients don't have to.



Mask Codes	
Code	Description
SLN1SMLZ	F&P Solo Nasal Mask Fit Pack - Small, Medium, Large Z Model
SLP1SMLZ	F&P Solo Pillows Mask Fit Pack - Small, Medium, Large Z Model
SLN1SZ	F&P Solo Nasal Mask Small
SLN1MZ	F&P Solo Nasal Mask Medium
SLN1LZ	F&P Solo Nasal Mask Large
SLN1WZ	F&P Solo Nasal Mask Wide
SLP1SZ	F&P Solo Pillows Mask Small
SLP1MZ	F&P Solo Pillows Mask Medium
SLP1LZ	F&P Solo Pillows Mask Large

Spare Parts	
Code	Description
400SL121	Solo Headgear
400SL141	Solo Tube, Frame and Swivel
400SLN111	Solo Nasal Cushion - Small
400SLN112	Solo Nasal Cushion - Medium
400SLN113	Solo Nasal Cushion - Large
400SLN114	Solo Nasal Cushion - Wide
400SLP111	Solo Pillows Cushion - Small
400SLP112	Solo Pillows Cushion - Medium
400SLP113	Solo Pillows Cushion - Large

References: 1. 85 out of 91 participants rated their experience with the F&P Solo Nasal or F&P Solo Pillows masks as simple or very simple. Internal validation study conducted with 91 participants in USA in 2022.
2. 79 out of 83 participants were able to set themselves up on the F&P Solo Nasal or F&P Solo Pillows masks without assistance. Defined as orientating the mask on the head/face as intended and achieving an acceptable seal without human interaction. Internal self-setup trial conducted with 83 participants, naïve to CPAP therapy, in New Zealand in 2022. F&P, AutoFit, AutoLock, F&P Solo and myMask are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip. 629105 REV B © 2024 Fisher & Paykel Healthcare Limited.