



Steam Challenge Pack is designed for steam sterilization procedures. The package contains Type 5 (Class 5) chemical integrator **(which gives equivalent results from biological indicators)** and high quality barrier papers for to measure penetration.

The front side of integrator has a window which contains **REJECT** and **ACCEPT** statements. If the acceptance line has been passed, it means many parameters of sterilizer have been tested and the test is considered successful.

- ★ It is easy to use.
- ★ Water-based and non-toxic chemical indicator is used.
- ★ It is easy to read and interpret the integrator result.
- ★ The test page has fields to fill in and can be kept as a permanent record for traceability. It's suitable for archiving.
- ★ Compatible with ISO 11140-4.

CP-200



STERISTAR Record Card Complies with
Type 5 Integrator Challenge Pack **STEAM** EN-450 11140-1
For use at 134°C, 3.5min. or 120°C for 9 min. or 121°C, 16.5min. TYPE 5 Integrator

Date: _____ Dept: _____
Sterilizer N°: _____ Load N°: _____
Temp: _____ Exposure Time: _____
Operator: _____ Approved or Not: _____

STERISTAR
TYPE 5
Complies with
EN4501140-1
STEAM REF CP-200
LOT 2010983
2020-10
2025-10

Dark Star Must Pass This Point

REJECT **ACCEPT**

Safe, non-toxic chemicals; no risk for contamination or transfer to sterile items.

REF CP-200
LOT 2010983

2020-10
2025-10

UNPROCESSED

CARD SIZE:
124 X 124 mm



UNPROCESSED



REJECTED



ACCEPTED

REF NO	DESCRIPTION	PACKAGING
CP-200	Steam Challenge Pack	80 Pcs. / CARTON