




Sterilization container labels are designed to be used in the pocket of the containers. Two kinds of container labels are available as self-adhesive and non-self-adhesive label. It is easy to understand whether sterilization has been performed because of the indicators on the label.


- ★ Water-based and non-toxic chemical indicator is used.
- ★ It is easy to read and interpret.
- ★ Compatible with ISO 11140-1.
- ★ Different designs are made for the company/person on request.
- ★ Self-Adhesive versions can be used as a permanent record for traceability purposes.

C-1075




PRODUCT SIZE: 67 X 37 mm

C-1075-SA



PRODUCT SIZE: 67 X 37 mm

C-1076



PRODUCT SIZE: 58 X 17 mm

REF NO	DESCRIPTION	PACKAGING
C-1075	Sterilization Container Label with Steam Indicator (Pack of 100)	60 Packs / CARTON
C-1075-SA	Self Adhesive Sterilization Container Label with Steam Indicator (Pack of 100)	60 Packs / CARTON
C-1076	Sterilization Container Label with Steam Indicator, Small (Pack of 100)	80 Packs / CARTON