


Sterilization Container Self Adhesive Seals



Container seal labels have been specially designed to safeguard the container from being opened by unauthorized personnel. The function of the container label is to safeguard and mark containers.

- ★ Water-based and non-toxic chemical indicator is used.
- ★ It is easy to read and interpret.
- ★ Compatible with ISO 11140-1.
- ★ Different designs are made for the company/person on request.
- ★ Self-Adhesive versions can be used as a permanent record for traceability purposes.
- ★ Non-Self-Adhesive versions are available on request.



PRODUCT SIZE: 140 X 30 mm



Unglue the tail section on the self-adhesive label.



Pass it through the lock section of the container.



The label is locked with the adhesive feature. For archiving, fill in the corresponding locations on the label.

REF NO	DESCRIPTION	PACKAGING
C-1074	Self Adhesive Container Seal with Steam Indicator (Pack of 100)	50 Packs / CARTON

