




A reliable sterilization tape that provides immediate identification of processed items and also serves as a closure tape to seal sterilization packs. Chemical indicator strips change color when exposed to appropriate sterilization conditions.

The sterilization tapes are suitable for steam, ethylene oxide and hydrogen peroxide sterilization processes.

- ★ Water-based and non-toxic chemical indicator is used.
- ★ It adheres to woven, treated woven, non-woven, paper and paper/plastic wraps.
- ★ The adhesive feature ensures that the winding material is safely packed and removed from the material in a clean/easy way.
- ★ Compatible with ISO 11140-1.

Sterilization tape changes color when subjected to sterilization process. It is important to note that the presence of sterilization tape that has changed color on an item does not ensure that the product is sterile, the tape only shows if sterilization has been performed.



| REF NO | DESCRIPTION | PACKAGING |
|----------|-----------------------------------|------------------|
| SAT-1950 | Steam Tape 19mm x 50M | 60 Pcs. / CARTON |
| SAT-2550 | Steam Tape 25mm x 50M | 50 Pcs. / CARTON |
| B-1950 | Steam Tape 19mm x 50M | 60 Pcs. / CARTON |
| B-2550 | Steam Tape 25mm x 50M | 50 Pcs. / CARTON |
| E-1950 | Ethylene Oxide Tape 19mm x 50M | 60 Pcs. / CARTON |
| E-2550 | Ethylene Oxide Tape 25mm x 50M | 50 Pcs. / CARTON |
| H-1950 | Hydrogen Peroxide Tape 19mm x 50M | 60 Pcs. / CARTON |
| NPT-1950 | Non-printed Tape 19mm x 50M | 60 Pcs. / CARTON |
| NPT-2550 | Non-printed Tape 25mm x 50M | 50 Pcs. / CARTON |