

# Hausted® Powered and Hydraulic Mobile Surgi-Stretcher for Ophthalmology



Hausted®

## Stretcherology: 5E8/578 Surgi-Stretcher Series

The Surgi-Stretcher Series is designed as a mobile surgical platform for EYE/ENT and plastic surgery procedures.

### Versatility

- Articulating head piece with dual operating knobs allow infinite adjustability and precise microscopic movement of the head section
- Designed with a low-profile head piece allowing site access for both superior (over brow) and lateral (side approach) procedures
- Easy-to-use electric/battery powered or hydraulic/foot-end operated controlled backrest and leg flex operated controls are available for height adjustment, backrest, leg section, and Trendelenburg/Reverse Trendelenburg



Microscopic articulation of the head piece is achieved by rotating the control knobs to raise/lower the patient's forehead or chin.

### Product Specifications

MODEL	ELECTRIC/BATTERY POWERED 5E8-EYE-ST 5D8-EYE-ST* 5E8-EYX-ST 5D8-EYX-ST*	HYDRAULIC/ FOOT-END CONTROLLED BACKREST AND LEG FLEX 578-EYE-ST 578-EYX-ST
Height Range: High	30 in [76 cm]	32 in [81 cm]
Height Range: Low	22 in [56 cm]	24 in [61 cm]
Overall Width	29.75 in [75.6 cm]	29.75 in [75.6 cm]
Overall Width (EYX Models Only)	34.5 in [87.6 cm]	34.5 in [87.6 cm]
Patient Surface Width	24.5 in [62.2 cm]	24.5 in [62.2 cm]
Patient Surface Width (EYX Models Only)	29.25 in [74.3 cm]	29.25 in [74.3 cm]
Overall Length	88.5 in [225 cm]	88.5 in [225 cm]
Patient Surface Length	81.5 in [207 cm]	81.5 in [207 cm]
Backrest Degree of Movement	0-80°	0-90°
Knee Flex Degree of Movement	0-60°	0-75°
Standard Trendelenburg	10°	16°
Retracto Rail Height	14 in [35.5 cm]	14 in [35.5 cm]
Retracto Rail Length	49 in [125 cm]	49 in [125 cm]
Maximum Weight Capacity Evenly Distributed	720 lbs [327 kg] * DC/220 volt version	500 lbs [227 kg]

### Standard Features

- 8 in [20 cm] easy rolling casters
- Four wheel brake and steer caster system operable from all four corners
- Easy to clean thermoformed base cover with O<sub>2</sub> tank holder and ample storage space
- Low profile articulating head piece
- Dual operating control knobs
- Stainless steel Retracto® Side Rails with rail caps
- 3 in [7.6 cm] comfort mattress and head pad
- 3-Stage Telescoping I.V. pole [permanently mounted at foot]
- Dual side hydraulic foot pedal controls [578]
  - Height Adjustment
  - Trendelenburg Positioning
- Foot-end hand crank controlled backrest and leg flex [578]

### Standard Features for Electric/Battery Powered [5E8-EYE-ST]

- Power hand held controls:
  - Height Adjustment
  - Trendelenburg/Reverse Trendelenburg
  - Backrest Adjustment
  - Knee Flex Adjustment
  - Auto Contouring Adjustment
- Nurse Control Panel at foot-end
- Emergency Backrest Release lever to override electric/battery
- Battery back-up allows power functions when not plugged in



770.368.4700

www.hausted.com

GF1500008RevB15



**GRAHAM-FIELD®**  
Proudly made in the USA

© 2015 GF Health Products, Inc